



Microloops® Loops and needles for perfect inoculation.

Microloops[®]


A comprehensive range of inoculating loops. Decades of experience in the manufacture of precision inoculating loops ensure you can rely on Microloops® for accurate volumetric technique every time.

Microloops® Plastic

Premium quality sterile plastic inoculating loops and needles. The precision moulded smooth edges prevent tearing of the agar surface, while good wettability ensures safe transfer of liquid samples. The Classic style features a dual sectional shaft, rigid but flexible for comfortable handling and control.

Precision moulded plastic loops and needles are available in individual peel pouches, or in convenient tamper evident resealable packs of 20.

- 3 loop sizes: ΙμΙ, 5μΙ, 10μΙ.
- Needle version.
- Colour coded for ready selection
- Ultra smooth no gouging
- Good wettability
- Dual profile shaft (on Classic style) for controlled handling square section for control of streaking direction round section for easy rotation
- Rigid, but not too rigid good for sputum
- Flexible, but not too flexible for controllable streaking
- Tamper-evident bags for guaranteed sterility



Microloops® Nichrome 5

A range of premium quality, accurately sized, twisted Nichrome 5 wire loops, available in 6 volumetric sizes. Nichrome 5 wire offers all the benefits of platinum, without the cost. With finely twisted shank for flexibility, the ultrasmooth precision formed loops glide smoothly across the agar surface without tearing, producing properly isolated colonies.

- Fully closed loops for secure transfer of drop
- Twisted wire reduced vibration
- Uniform length
- Nichrome 5 stable during repeated heat/cool cycles
- Reduced oxidation
- Reduced carbonisation



C.H. Collins, MBE FIMLS MIBiol, 1976

"The Medical Wire & Equipment Company loops are about the right length and, being twisted, do not vibrate as much as single wires. We have found them entirely satisfactory." \(^1\)

Standard or Calibrated

Standard Nichrome 5 Microloops® are packed in 25's, are carefully and accurately manufactured, and are completely suitable for consistent and precise transfer of sample materials. Each batch is carefully checked to ensure the loops meet the specification for the particular product.

Calibrated Nichrome 5 or Platinum Microloops® are individually packed, and each loop has been individually measured to ensure it meets the specification. A calibration certificate is supplied with each loop.

Microloops® made from Nichrome 5





Nichrome 5 is an alloy of 80% Nickel and 20% Chromium and has long been regarded as the finest alloy of its type for repeated

heating. It gives outstanding performance at temperatures up to 1200°C.

The Platinum Standard



The method streaking out of bacteria on agar surfaces plates was devised in Robert Koch's laboratory, but the introduction of platinum wire loops is attributed

to Mazyck P. Ravenel. ^{2,3} Although now largely superceded by nichrome, platinum is still used in some laboratories because of its faster cooling and complete resistance to oxidation and chemical degradation.

For those special applications use Medical Wire's Calibrated Platinum Microloops[®] available in 1 µl and 10 µl sizes. A straight wire version is also available.



platinum

A precious silvery-white metal, the chemical element of atomic number 78. Used in jewellery, electrical contacts, laboratory equipment, and industrial catalysts (Symbol: Pt).

Context & Method of Use

While the use of wire inoculating loops and needles dates back to the earliest days of bacteriology, the use of loops as measuring devices began in the dairy industry,⁵ and it was in the early 1960's that the inoculating loop began to be used as an actual diagnostic tool. O'Sullivan⁶ first described the culture of a "standard loopful" of urine, and McGeachie & Kennedy⁷ modified the technique to the now familiar streaking in successive planes.

The method of use has generated some debate over the years. The most current is to dip the end of the loop to just below the surface of the liquid (e.g. urine), and remove vertically, taking care not to carry over any liquid on the shank.⁸



Order Information

	٥	
-	<u></u>	?
-	_	•
_	0	-
	۵	ķ
	_	-

Microloops® P	Microloops® Precision Moulded Plastic						
Cat No.	Size	Style	Pack Size				
MW710/20	10μΙ	Classic	50 × 20	0-			
MW711/20	lμl	Classic	50 × 20	0			
MW705/20	5µl	Standard	50 × 20	0			
MW702/20	Needle	Standard	50 × 20				

Nichrome 5

Microloops® F	Microloops® Precision Nichrome 5 (Reusable)						
Cat No.	Size	Internal Diameter	Pack Size				
MW190	10µl	5.0mm	25	0			
MWI9I	5µl	3.5mm	25	0			
MW192	3.3µl	2.9mm	25	C			
MW193	2.5µl	2.5mm	25	0-			
MW194	2µl	2.2mm	25	0-			
MW195	lμl	1.5mm	25	D			
MW197	Straight W	/ire -	25	-			

Nichrome 5

Microloops® Calibrated Nichrome 5 (Reusable)		Calibrated		
Cat No.	Size	Internal Diameter	Pack Size	
MW1901	Ι0μΙ	Ø 5.0mm	1	0
MW1951	Lul	Ø I 5mm		



Microloops® Calibrated Platinum (Reusable)		Calibrated		
Cat No.	Size	Internal Diameter	Pack Size	
MW190PL	10μΙ	Ø 5.0mm	1	0
MW195PL	lμl	Ø 1.5mm		C-
MW197PL	Straight W	/ire	T	

Microloops® Holders

Ś
≍
<u></u>
C

Cat No.	Size	Pack Size	
MW196/6	150mm (6 Inch)	1	
MW196/7	175mm (7 Inch)	1	
MW196/8	205mm (8 Inch)	1	



Microloops [®]	Calibration Gauge	Calibrated	
Cat No.	Size	Pack Size	
MW188	Ι0μΙ		
MW189	lμl		

Calibration Calibration of Microloops® and Microstreakers® is carried out by the Evans Blue dye method and the FDA Drill Bit Method. 9

Precision manufacturing The manufacturing process ensures a perfectly symmetrical loop with equal dimensions, not an uneven oval. This allows the use of our gauge for future calibration once we have certified the loop.







Microstreakers[®]

Combined twisted Nichrome 5 wire loops and light weight, insulated aluminium holders, perfectly balanced for easy manipulation. Microstreakers glide smoothly across the surface of agar without tearing or gouging.

- Colour coded
- Flexible
- Non-gouging
- Rapid Cooling
- Consistent and reproducible results

Microloops® Calibration Gauge

Calibrated loops in regular use should be checked on a monthly basis for accuracy. An easy way to comply with this requirement is to use the Microloops® Calibration Gauge (UKAS certified). This is a simple, but precision tooled Go No-Go Gauge which is used in accordance with the US FDA drill bit method®. Measurement is physical and there is no requirement for a spectrophotometer.

- Go No-Go gauge for guaranteed calibration
- Simple to use
- IµI and I0µI
- Robust construction
- FDA drill bit method⁹



L-shaped Spreader

- Precision moulded plastic
- 30mm contact edge
- Convenient inoculation of larger volumes
- High visibility bright orange

CE-marking

Medical Wire's Microloops® products are CE-marked and meet the requirements for In Vitro Diagnostic Medical Devices as defined in the European In Vitro Diagnostic Devices Directive 98/79/EC. All Microloops® products are manufactured to ISO 9001:2008

Order Information

	Microstreakers® Precision Nichrome 5 (Reusable)					
	Cat No.	Size	Pack Size			
	MW180	Large	5			
Sio	MW182	Medium	5			
Precision	MW184	Small	5			
₫	MW186	Needle	5	===		

	Microstreakers® Calibrated Nichrome 5 (Reusable) Cat No. Size Pack Size MW1851 10µ1 1				
	Cat No.	Size	Pack Siz	ze	
dil	MW1851	Ι0μΙ		0	-
Ü	MW1871	lμl		L-	

Micro	loops [®] Precision I	Moulded Plastic - Indi	ividually Wrapped
Cat N	o. Size	Pack Size	
MW70	00/I I0µI	1000 × 1	0
MW70	01/I IµI	1000 × 1	0
MW70	05/I 5µI	1000 × 1	0
MW70	02/I Needle	1000 × 1	

er	L Shaped S	Spreader	
eac	Cat No.	Size	Pack Size
र्ट्	MW7121	30mm contact edge	500 × I

References

- 1. Collins, C.H., 1976 The bacteriologist's loop: a biohazard, The Gazette, 20:55-6
- Poupard, J.A., 2010, A History of Microbiology in Philadelphia: 1880 2010, Xlibris Corporation
- Ravenel, M. P., & J. Willoughby Irwin, 1905, Studies In Mixed Infection In Tuberculosis: A Preliminary Report in National Association For The Study And Prevention Of Tuberculosis Transactions Of The First Annual Meeting Washington, D. C. May 18th And 19th, 1905
- 4. Oxford Dictionaries. www.oxforddictionaries.com
- Burri, P, 1928, The quantitative smear culture: a simple means for the bacteriological examination of milk. Report of Proceedings World's Dairy Congress p690-697
- 6. O'Sullivan, D.J., et al, 1960, A simplified method for the quantitative bacterial culture of urine. J. Clin. Path., 13:527-8
- McGeachie, J., & A.C. Kennedy, 1963. Simplified quantitative methods for bacteriuria and pyuria. J. Clin. Path., 16:32-38
- Quality Guidance Q5, 2013, Inoculation of Culture Media for Bacteriology. UK Standards for Microbiology Investigations, Standards Unit, Microbiology Services, Public Health England, www.hpa.org.uk/SMI
- Isenberg H. (Editor in Chief) 2004 Clinical Microbiology Procedures Handbook -2nd edition.

