

## Insert for ESBL + AmpC Screen Kit (98008)

REVISION:	DBV0032G
DATE OF ISSUE:	15/03-2013
LANGUAGE:	English

### ESBL + AmpC Screen Kit

FOR IN VITRO DIAGNOSTIC USE ONLY

**PRODUCT GROUP:** Kits for beta lactamase identification

**MANUFACTURER:** ROSCO Diagnostica A/S, Taastrupgaardsvej 30, DK-2630 Taastrup, Denmark.

**INTENDED USE:** Tablets are used for *in vitro* identification of microbial resistance mechanisms by the agar tablet/disc diffusion method, in order to confirm the mechanism by which the organism has gained resistance to specific antimicrobial agents.

**INTENDED USERS:** Only to be used by professionals and people trained to work with microbes and disc diffusion testing.

#### **DETAILED**

#### **INSTRUCTIONS:**

ROSCO's detailed Instruction for Detection of Resistance Mechanisms using Neo-Sensitabs™ and Diatabs™ should be available in each laboratory working with ROSCO's products.

Latest edition of this book can be seen in and/or printed out from ROSCO's website [www.rosco.dk](http://www.rosco.dk).

Instructions for Use and User's Guide can be obtained free of charge from your local distributor on request, or from ROSCO Diagnostica A/S:

E-mail: [info@rosco.dk](mailto:info@rosco.dk) or

Fax: +45 43 52 73 74

#### **PRINCIPLE OF THE TEST:**

This combination screen test is designed for the routine screening of ESBL and/or AmpC producers. ESBL + AmpC Screen Kit consists of four cartridges of disc diffusion tablets: one cartridge of tablets with Cefotaxime alone and three cartridges of tablets with Cefotaxime combined with either Clavulanate (ESBL inhibitor) or Cloxacillin (AmpC inhibitor) and Cefotaxime combined with both Clavulanate and Cloxacillin. If an organism is suspected of ESBL and/or AmpC activity it can be shown by a difference in the inhibition zone of the cephalosporin alone and in combination with the inhibitors. The choice of cephalosporin is made to have the highest sensitivity and specificity with a single agent.

For the highest sensitivity and specificity the laboratory should include results for all third generation cephalosporins. All isolates with reduced susceptibility to these, or Aztreonam should be suspected of being ESBL producers.

For confirmatory AmpC and ESBL tests ROSCO recommends the use of our kits intended for these purposes; AmpC Confirm Kit and Total ESBL Confirm Kit.

**CONTENT AND FORMULATION:**

4 cartridges, each containing approximately 50 tablets, formulated for maximum stability:

- Cefotaxime 30 µg, coded CTX30
- Cefotaxime 30 µg + Clavulanate, coded CTX+C
- Cefotaxime 30 µg + Cloxacillin, coded CTXCX
- Cefotaxime 30 µg + Clavulanate + Cloxacillin, coded CTXCC

**STORAGE/HANDLING:**

Store at 2-8°C in the box provided or unopened cartridges until the expiry date shown on the product label. Allow the cartridges to acclimatize to room temperature for 30-60 minutes before the lid is removed from the cartridge. Once a cartridge has been opened and in particular when placed in a dispenser, it should be kept at room temperature for up to 2 months. If necessary, when in use for a longer period than 2 months, the cartridges can be stored at 2-8°C. Always seal the cartridges with the original green lid, and never place the dispenser in the refrigerator. When stored at 2-8°C the cartridges should be allowed to acclimatize, as described above, before use.

**PRECAUTIONS:**

For *in vitro* diagnostic use only. Safety precautions should be taken and aseptic techniques used when working with potential biohazards. To be used only by adequately trained and qualified laboratory personnel. Sterilize all biohazard waste before disposal. Refer to Product Safety Data Sheet.

**MATERIALS REQUIRED BUT NOT PROVIDED:**

Standard microbial equipment such as loops, culture media, incubator etc. and biochemical reagents.

**PROCEDURE:**

1. Using a fresh, pure culture prepare a suspension of the organism to be tested equivalent to McFarland 0.5
2. Using a sterile swap or Drigalski spatula spread the suspension uniformly over the entire area of a Mueller Hinton susceptibility agar plate
3. Using a single tablet dispenser, place one of each tablet on the inoculated agar plate, ensuring sufficient space between individual tablets to allow for proper measurement of inhibition zones. Notice that more than one Screen Kit can be tested on the same plate.
4. Incubate at 35±1°C for 18±2 hours (overnight)
5. Measure and record the diameter of the inhibition zone. No zone around a tablet corresponds to a 9 mm inhibition zone.

**INTERPRETATION OF RESULTS:**

The results are interpreted by comparing the inhibition zones of the different tablets.

1. Compare the zone of inhibition of the Cefotaxime 30 µg tablet to the zones of inhibition of each of the Cefotaxime 30 µg + inhibitor(s) combination tablets. If all zones are within 3mm of each other, record the organism as neither expressing ESβL nor AmpC activity.
2. Measure the inhibition zones around Cefotaxime 30 µg (CTX30) and Cefotaxime 30 µg + Cloxacillin(CTXCX) and compare with Cefotaxime 30 µg + Clavulanate(CTX+C) and Cefotaxime 30 µg + Clavulanate + Cloxacillin (CTXCC), respectively. If both CTX+C – CTX30 and CTXCC – CTXCX is ≥ 5mm and both the values in step 3) are <5 mm, the organism is demonstrating ESβL activity alone.
3. Measure the inhibition zones around Cefotaxime 30 µg + Clavulanate (CTX+C) and Cefotaxime 30 µg (CTX30) and compare with Cefotaxime 30 µg + Clavulanate + Cloxacillin(CTXCC) and Cefotaxime 30 µg + Cloxacillin(CTXCX), respectively. If both CTXCC – CTX+C and CTXCX – CTX30 ≥ 5mm and the values in step 2) is < 5mm, the organism is demonstrating AmpC activity alone.
4. Measure the inhibition zones around Cefotaxime 30 µg (CTX30) and Cefotaxime 30 µg + Cloxacillin (CTXCX) and compare with Cefotaxime 30 µg + Clavulanate(CTX+C) and Cefotaxime 30 µg + Clavulanate + Cloxacillin(CTXCC), respectively. If CTXCC – CTXCX ≥ 5mm (ESBL) and CTX+C – CTX30 < 5mm (AmpC), the organism is demonstrating both ESβL and AmpC activity.  
Please notice that Clavulanate has two different functions:
  1. It inhibits ESBL and
  2. It induces AmpC. This explains why CTX+C-CTX30 < 5mm, if the isolate possesses an AmpC
5. Use table 1 to assist in the interpretation

**Procedure for reading the results.**

- Step 1:** If B-A and/or D-C ≥ 5 mm: ESBL positive.  
Otherwise ESBL is negative.
- Step 2A:** Isolate found ESBL positive:  
If D-B ≥ 5 mm and/or B-A < 5 mm: AmpC positive.  
Consequently isolate ESBL + AmpC positive
- Step 2B:** Isolate found ESBL negative:  
If C-A and/or D-B >= 5 mm: AmpC positive.  
Consequently isolate ESBL negative and AmpC positive  
Otherwise negative for both ESBL and AmpC

A calculator program is available.

**QUALITY CONTROL:**

Although ROSCO Diagnostica A/S produces the, by far, most stable diffusion discs (tablets) it is necessary to perform regular quality control. This should be done with at least one organism to demonstrate a positive reaction and at least one organism

to demonstrate a negative reaction. Zones of inhibition obtained using the combination tablets plus the cephalosporin alone tablet against the negative control (i.e. *E. coli* ATCC 25922), should be within 3 mm. Any greater difference indicates that the product has lost activity and should not be used.

As positive Q. C strains the following may be used:  
*Enterobacter cloacae* NCTC 13406, Amp C positive  
*Enterobacter cloacae* ATCC BAA – 1143, Amp C positive  
*Klebs. pneumoniae* ATCC 700603 ESBL positive

**Table 1**

		<b>Cefotaxime CTX30 A</b>	<b>Cefotaxime+Clav. CTX+C B</b>	<b>Cefotaxime+Cloxa. CTXCX C</b>
<b>ESβL</b>	CTX+C or B CTXCC D	≥ 5 mm -	- <5 mm	- ≥ 5 mm
<b>AmpC</b>	CTXCX or C CTXCC D	≥ 5 mm -	- ≥ 5 mm	- <5 mm
<b>ESβL + AmpC</b>	CTX+C and B CTXCC D	<5 mm (AmpC) -	- ≥ 5 mm (AmpC)	- ≥ 5 mm (ESBL)

Neither ESβL nor AmpC: All zones within 3mm of each other. CTXCC = Cefotaxime + Clavulanate + Cloxacillin

Note : “-“ means that the difference is irrelevant for the mechanism (i.e. the difference between CTXCC compared to CTX30 is irrelevant for the detection of ESBL).

REFERENCES: [www.rosco.dk](http://www.rosco.dk)