

UREASE TEST TABLETS

DISCUSSION:

Several species of Enterobacteriaceae, including *Proteus*, *Klebsiella*, and some *Citrobacter* species, as well as some *Haemophilus* species and several other bacteria, fungi, and yeast species produce the enzyme urease. This enzyme hydrolyzes urea into ammonia, water and carbon dioxide. (1) Since most stool pathogens, particularly *Salmonella*, *Shigella*, and *E. coli*, are urease negative, it can be useful to determine an organism's urease activity. A notable exception to this is *Yersinia enterocolitica*, which is urease positive. Key Urease Test tablets indicate the presence of the alkaline ammonia end product by a color change from straw color to pink or red.

MATERIAL SAFETY:

Each K650 tablet contains 10 mg. of urea. K1650

tablets contain less than 2 mg. of urea. This product does not contain any material known at this time to be hazardous.

MATERIALS:

K650 Urease Tablets are sold on bottles of 50 tablets only. K1650 WEE-TABS are sold in packs of 28 ready-to-use, pre-labeled tubes. Urease test tablets require preliminary growth on media appropriate for the specimen. Consult a reference manual such as the Manual of Clinical Microbiology for the appropriate media to use. The following items are required but not provided:

Test tubes (K650 only)

Loop for harvesting colonies

Distilled water, neutral pH

QUALITY CONTROL:

Key Urease tablets should be tested prior to use with organism of known reactivity. Key suggests *Proteus vulgaris*

ATCC 8427 as positive and *E. coli* ATCC 25922 as negative. Discard used tests in a manner conforming with accepted laboratory procedure for bio-hazardous materials.

INSTRUCTIONS:

1) Add 1 Urease Test Tablet to 1 ml. of distilled water in a small test tube. For WEE-TAB Urease add 3-5 drops of water.

2) Inoculate heavily from a culture grown on a slant or plate. Kligler, Russell, or T.S.I. agar may be used without affecting the results.

3) Incubate at 35-37C. *Proteus* spp give a positive test in 1-8 hours. Some *Klebsiella* strains produce urease in smaller amounts and may require up to 24 hours in the K650 to show a positive reaction. Any pink or red color is a positive result.

WEE-TAB reactions take as little as 15 minutes and are finished at 6 hours.

REACTIONS:

Organism	Reaction	% Positive
<i>Citrobacter</i> sp.	+	78
<i>Edwardsiella</i>	-	0
<i>Enterobacter</i> sp. *	-	10
<i>E. coli</i>	-	1
<i>Hafnia</i>	-	0
<i>Kleb. pneumoniae</i>	+	96
<i>Morganella</i>	+	98
<i>Proteus mirabilis</i>	+	98
<i>Proteus vulgaris</i>	+	98
<i>Providencia</i> sp.*	+/-	30
<i>Salmonella</i> sp.	-	99
<i>Shigella</i> sp.	-	0
<i>Yersinia</i> sp. *	+	80

* Reactions vary within species

REFERENCES:

(1) Bailey and Scott's Diagnostic Microbiology, 8th Editions, Chapters 9 and 27.



KEY SCIENTIFIC PRODUCTS, INC
1113 EAST REYNOLDS STREET
STAMFORD, TEXAS 79553
VOICE 800-843-1539
FAX 888-440-4208
WWW.KEYSCIENTIFIC.COM