



KEY SCIENTIFIC PRODUCTS
1113 EAST REYNOLDS ST.
STAMFORD, TEXAS 79553
WWW.KEYSCIENTIFIC.COM
Voice 800-843-1539
Fax 888-440-4208

Safety Data sheet

1 Identification

GHS Product Identifier

Catalog Number	Product Name: Multiple
K1048	WEE-TAB PO4/MAGL/ARG
K1083	WEE-TAB Arginine Naphthylamide
K1084	WEE-TAB Trypsin Naphthylamide
K1086	WEE-TAB Chymotrypsin Naphthylamide
K1135	WEE-TAB BGLU/PYR
K1137	WEE-TAB ESC/PYR
K1139	WEE-TAB Strep Bovis Screen
K1202	WEE-TAB CEAL (w/PEP Reagent)
K1236	WEE-TAB AFU/LGY
K1237	WEE-TAB AFU/MAGA/CHY
K1242	WEE-TAB AGLU/PAL
K1272	WEE-TAB AGAL/SER
K1276	WEE-TAB BGLU/ARG
K1277	WEE-TAB AGLU/ARG
K1305	WEE-TAB Glycine Naphthylamide
K1378B	LAP Discs (w/PEP Reagent)
K1462	WEE-TAB NAG/Proline Naphthylamide
K1466	WEE-TAB NAG/MBGA/TRY
K1495	WEE-TAB ONPG/PYR
K1525	WEE-TAB Phenylalanine Naphthylamide
K1532B	PRO Discs (w/PEP Reagent)
K1538B	PYR Discs (w/PEP Reagent)
K1538N	PYR Strips (w/PEP Reagent)
K1539	COLI ID Discs
K880481	Nitrate B 30 mL Dropper Bottle
K9140	WEE-TAB NGP
K9142	Candida Screen (w/Reagents)
K9131B	Butyrate Discs (w/PEP Reagent)
K9162B	Try Discs (w/PEP Reagent)
K9187	WEE-TAB Neisseria
K9265B	GGT Discs (w/PEP Reagent)
K971538	SWABIT PYR
K9760	Clostridium Screen Discs
K979760	SWABIT Clostridium Screen
K980481	DROPIT Nitrate B
KF4012	Identicult AE Cards (w/Reagents)
KF4020	Identicult PRP Card (w/Reagents)
KF4080	Identicult AE Strips (w/Reagents)
KS5300	Fluo-Card Milleri
KS5400	Flow-Card Strep

Other means of identification

K1048	P04/MAGL/ARG
K1083	Arginine

K1084	Trypsin
K1086	Chymotrypsin
K1135	BGLU/PYR
K1137	ESC/PYR
K1139	Strep Bovis
K1202	Ceal
K1236	AFU/LGY
K1237	AFU/MAGA/CHY
K1242	AGLU/PAL
K1272	AGAL/SER
K1276	BGLU/ARG
K1277	AGLU/ARG
K1305	Glycine
K1378B	LAP Discs
K1462	NAG/Proline
K1466	NAG/MBGA/TRY
K1495	ONPG/PYR
K1525	Phenylalanine
K1532B	Proline
K1538B	PYR Discs
K1538N	Pyr Strips
K1369	Coli ID
K880481	Nitrate B Dropper
K9140	NGP
K9142	Candida
K9131B	Butyrate Discs
K9162B	Try Discs
K9187	Neisseria
K9265B	GGT Discs
K971538	PYR Swab
K9760	Clostridium Discs
K979760	Clostridium Swabs
K980481	Nitrate B Droppers
KF4012	Identicult AE Cards
KF4020	Identicult PRO Cards
KF4080	Identicult AE Strips
KS5300	Fluo-Card Milleri
KS5400	Flow-Card Strep

Recommended use of the chemical and restriction on use

For invitro diagnostic use only by trained professionals.

Supplier's details

Manufacturer / Supplier:

Key Scientific Products, Inc.
1113 East Reynolds Street
Stamford, TX 79553

Phone Number: 1-800-843-1539
Emergency Phone Number: none available.

2 Hazard(s) identification

Classification of the substance or mixture

Carcinogenicity (Category 1B,2), H350, H351
Acute toxicity, Oral (Category 4),H302

Skin irritation (Category 2), H315

Eye irritation (Category 2A), H319

Specific organ toxicity - single exposure (Category 3), Respiratory system, H335

Acute aquatic toxicity (Category 1,2), H400, H401

Chronic aquatic toxicity (Category 1, 2), H410, H411

GHS label elements

Warning



Toxic if swallowed

Harmful if swallowed

Causes skin irritation

Causes serious eye irritation

May cause respiratory irritation

Suspected of causing cancer

Toxic to aquatic life

Toxic to aquatic life with long lasting effects

Obtain special instructions before use.

Do not handle until all safety precautions have been read and understood.

Avoid breathing dust/fume/gas/mist/vapours/spray.

Wash skin thoroughly after handling.

Do not eat, drink or smoke when using this product.

Use only outdoors or in a well-ventilated area.

Avoid release to the environment.

Wear protective gloves/protective clothing/eye protection/face protection.

Use personal protective equipment as required.

IF SWALLOWED: Call a POISON CENTER/ doctor/physician if you feel unwell.

IF ON SKIN: Wash with plenty of soap and water.

IF INHALED: Remove person to fresh air and keep comfortable for breathing.

IF exposed or concerned: Get medical advice/ attention.

If skin irritation occurs: Get medical advice/attention.

If eye irritation persists: Get medical advice/attention.

Wash with plenty of soap and water.

Take off contaminated clothing and wash before reuse.

Collect spillage.

Store according to label instructions.

Store locked up.

Dispose of contents/container to an approved waste disposal site.

Other hazards which do not result in classification

Follow directions on label and technical insert.

3 Composition/information on ingredients

Description	CAS Number	EINECS Number	%	Note
L-Glutamic acid γ -(β -naphthylamide)	14525-44-1	238-542-7	0.1 - 0.5	
L-Leucine β -naphthylamide hydrochloride	893-36-7	212-970-4	0.1 - 0.5	
L-Proline β -naphthylamide hydrochloride	97216-16-5		0.1 - 0.5	
L-Pyroglutamic acid 2-naphthylamide	22155-91-5	244-809-9	0.1 - 0.5	
Na-Benzoyl-DL-arginine β -naphthylamide hydrochloride	913-04-2	213-015-4	0.1 - 0.5	
L-Phenylalanine β -naphthylamide	740-57-8		0.1 - 0.5	
L-Serine β -naphthylamide	888-74-4		0.1 - 0.5	
Glycine β -naphthylamide	716-94-9		0.1 - 0.5	
4-Nitrophenyl N-acetyl- β -D-glucosaminide	3459-18-5	222-398-7	0.1 - 0.5	
2-Nitrophenyl β -D-galactopyranoside	369-07-3	206-716-1	0.1 - 0.5	
N,N-Dimethyl-1-naphthylamine	86-56-6	201-682-4	0.1 - 0.5	
N-Glutaryl-Gly-Gly-Phe- β -naphthylamide	55145-12-5		0.1 - 0.5	
trans-4-Hydroxy-L-proline β -naphthylamide	3326-64-5	222-045-7	0	

4 First-aid measures

Description of necessary first-aid measures

Eyes: In case of contact with eyes, rinse immediately with water for 10-15 minutes. If irritation develops and persists seek medical advice.

Ingestion: Drink large amounts of water. Seek medical attention.

Inhalation: Not applicable. No vapors.

Skin: Wash thoroughly with soap and water. If irritation occurs seek medical advice.

Most important symptoms/effects, acute and delayed

No data available.

Indication of immediate medical attention and special treatment needed, if necessary

No data available.

5 Fire-fighting measures

Suitable extinguishing media

Use waterspray, CO₂, foam, or dry powder as the extinguisher medium.

Specific hazards arising from the chemical

Carbon oxides, nitrogen oxides, hydrogen chloride gas.

Special protective actions for fire-fighters

Wear self-contained breathing apparatus if needed.

6 Accidental release measures

Personal precautions, protective equipment and emergency procedures

Wear lab coat, gloves, and safety glasses.

Methods and materials for containment and cleaning up

Sweep or vacuum uninoculated product and discard as normal waste.

Inoculated tests should be treated according to accepted procedures for bacteriological hazards.

7 Handling and storage

Precautions for safe handling

Wear appropriate personal protective equipment (gloves, lab coat, safety glasses).

Conditions for safe storage, including any incompatibilities

Store according to label directions in the sealed container provided.

8 Exposure controls/personal protection

Control parameters

Contains no substances with occupational exposure limit values.

Appropriate engineering controls

Handle in accordance with good industrial hygiene and safety practices. Wash hands before breaks and at the end of workday.

Individual protection measures

Eye/face protection

Safety glasses with side shields conforming to EN166. Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or en 166 (EU).

Skin protection

Handle with gloves. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands.

Body protection

Choose body protection in relation to its type, to the concentration and amount of dangerous substances, and to the specific work-space.

Respiratory protection

Respiratory protection is not required. Where protection from nuisance levels of dusts are desired, use type N95 (US) or type p1 (EN 143) dust masks. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (UN).

Control of environmental exposure

Do not let product enter drains.

9 Physical and chemical properties

Physical and chemical properties

Appearance:

Form: Solid

Color: White disc or white/gray tablet.

Odor:	No data available.
Odor Threshold:	No data available.
pH:	No data available.
Melting point/freezing point:	No data available.
Melting point:	No data available.
Flash Point:	No data available.
Evaporation rate:	No data available.
Flammability:	No data available.
Upper/lower flammability/explosion limits:	No data available.
Vapor Pressure:	No data available.
Vapor density:	No data available.
Relative density:	No data available.
Water solubility:	Soluble in water.
Partition coefficient (n-octal/water):	No data available.
Auto-ignition temperature:	No data available.
Decomposition temperature:	No data available.
Viscosity:	No data available.
Explosive Properties:	No data available.

10 Stability and reactivity

Reactivity

No data available.

Chemical stability

Stable under recommended storage conditions.

Possibility of hazardous reactions

No data available.

Conditions to avoid

No data available.

Incompatible materials

Oxidizing agents.

Hazardous decomposition products

No data available.

In the event of a fire, see Section 5.

11 Toxicological information

Toxicological (health) effects

Acute toxicity:

No data available.

Inhalation: No data available.

Dermal: No data available.

Skin corrosion/irritation

No data available.

Serious eye damage/eye irritation

No data available.

Respiratory or skin sensitisation

No data available.

Germ cell mutagenicity

No data available.

Carcinogenicity

IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible, or confirmed human carcinogen by IARC.

ACGIH: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible, or confirmed human carcinogen by ACGIH.

NTP: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible, or confirmed human carcinogen by NTP.

OSHA: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible, or confirmed human carcinogen by OSHA.

Reproductive toxicity

No data available.

Specific organ toxicity - single exposure

No data available.

Specific organ toxicity - repeated exposure

No data available.

Additional information

RTECS: Not available.

These chronic effects are suggested from analogy to beta-naphthylamine, RTECS# QM2100000. To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

12 Ecological information

Toxicity

No data available.

Persistence and degradability

No data available.

Bioaccumulative potential

No data available.

Mobility in soil

No data available.

Results of PBT and vPvB assessment

PBT/vPvB assessment not available as chemical safety assessment not required/not conducted.

Other adverse effects

No data available.

13 Disposal considerations

Disposal methods

Product

Uninoculated items may be discarded as normal waste.

Inoculated waste should be discarded in a manner appropriate for biological hazards.

Contaminated Packaging

Dispose of as unused product.

14 Transport information

UN Number

DOT (US)

Not dangerous goods.

IMDG

Not dangerous goods.

IATA

Not dangerous goods.

15 Regulatory information

Safety, health and environmental regulations specific for the product in question

SARA 302 Components

SARA 302: No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

SARA 313 Components

SARA 313: This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimus) reporting levels established by SARA Title III, Section 313.

SARA 311/312 Hazards

No SARA Hazards.

Massachusetts Right To Know Components

No components are subject to the Massachusetts Right to Know Act.

Pennsylvania Right To Know Components

No components are subject to the Pennsylvania Right to Know Act.

New Jersey Right To Know Components

No components are subject to the New Jersey Right to Know Act.

California Prop. 65 Components

These products do not contain any chemicals known to State of California to cause cancer, birth defects, or any other reproductive harm.

16 Other information

Other information

The above information, to the best of our knowledge, is accurate. Key Scientific Products assumes no liability whatsoever for the accuracy or completeness of the information stated above. Final determination of suitability of any material is the sole responsibility of the user. All materials may present unknown hazards and should be used with caution. Although certain hazards may be described, we cannot guarantee that these are the only hazards that exist.