

KEY SCIENTIFIC PRODUCTS

Serving Science Since 1935



MYCOBACTERIA DIGESTANT KITS AND COMPONENTS

Because specimens containing mixed flora of mycobacteria and other organisms makes recovery difficult, the use of NALC digestants/decontaminants is common practice in mycobacteriology laboratories worldwide.

Key Scientific's NAOH-NALC Digestant Kits and reagents are conveniently packaged and offer long shelf life.

In addition, NAOH-NALC RED's indicator makes it easy to determine the correct pH when adding phosphate buffer after digestion/decontamination.



CATALOG # - DESCRIPTION	UNIT	LIST
K7620050 NAOH-NALC KIT 2.0% 50 ML	5/BOX	\$34.00
K7620200 NAOH-NALC KIT 2.0% 200 ML	5/BOX	\$34.00
K7625050 NAOH-NALC KIT 2.5% 50 ML	5/BOX	\$34.00
K7625200 NAOH-NALC KIT 2.5% 200ML	5/BOX	\$110.00
K7640050 NAOH-NALC KIT 4.0% 50 ML	5/BOX	\$34.00
K7640200 NAOH-NALC KIT 4.0% 200ML	5/BOX	\$110.00
K7670050 PHOSPHATE BUFFER 50ML	25/BOX	\$87.00
K7670500 PHOSPHATE BUFFER 500 ML	EACH	\$19.00
K7670505 PHOSPHATE BUFFER 500 ML	10PK	\$159.00
K7680003 BSA (BOVINE SERUM ALBUMIN)3ML	20/BX	\$32.00
K7R25050 NAOH-NALC RED KIT 2.5% 50 ML	5/BOX	\$34.00
K7R25200 NAOH-NALC RED KIT 2.5% 200 ML	5/BOX	\$110.00

1113 E. Reynolds Street

Stamford, TX 79553

Phone: 800-843-1539

FAX: 888-440-4208

www.keyscientific.com