


Reagent for Measuring α -fetoprotein by latex-enhanced turbidimetric immunoassay

AFP-LATEX "SEIKEN"


Excellent sensitivity
Dynamic measuring range : 7 - 250 ng/mL

REFERENCE DATA (Hitachi 917)

Lower detection limit



Linearity




Within-run precision

	Low	Medium	High
n	20	20	20
Mean	9.9	22.6	99.5
SD	0.4	0.31	0.69
%CV	4.0%	1.4%	0.7%

REFERENCE DATA (Hitachi 917)

Prozone



Correlation with RIA

