

Reagent for Measuring Ferritin levels by latex-enhanced turbidimetric immunoassay

FER-LATEX X2 "SEIKEN" CN



Excellent performance for anemia and cancer testing

Dynamic measuring range: <4 - 1000 ng/mL



High sensitivity & Wide range

REFERENCE DATA (Hitachi 917)

Sensitivity



Linearity & Prozone



Precision

	Within-run			Between-day		
	Low	Mid	High	Low	Mid	High
N	21	21	21	20	20	20
MEAN	14.9	100.0	431.1	16.5	105.2	428.7
SD	0.60	0.65	2.20	0.87	1.60	3.52
%CV	4.0	0.6	0.5	5.3	1.5	0.8

