



### PLASTIC MICROLOOPS

**Description:** Premium Quality sterile plastic inoculating loops and needles. Precision moulded plastic microloops and needles are available in individual peel pouches, or in convenient tamper evident resealable packs of 20.

**Sizes available:** 1µl, 5µl and 10µl, plus needle version.

**Opening the bags:** Tear off the perforated section, initially try to gently pull the two layers apart to open the Ziplock, Note: if the package does not open easily, slide the Ziploc horizontally between two fingers to open the package.

**Method of use:** Remove loop from sterile plastic bag, and collect test material such as urine specimen, or bacterial colonies from an agar plate. Once used dispose of loop as clinical waste. Re seal the bags in the normal fashion by sliding fingers across the ziplock.

**Storage: No specific Storage Temperatures**

**Expiry Date: This product has a five year expiry date.**



Medical Wire & Equipment Co (bath) Ltd  
Corsham  
Wiltshire SN13 9RT UK  
[info@mwe.co.uk](mailto:info@mwe.co.uk) [www.mwe.co.uk](http://www.mwe.co.uk)



Revision: 4 Date 10/03/2020



Advena Ltd, Tower Business Centre, 2<sup>nd</sup> Fl.,  
Tower Street, Swatar, BKR 4013, Malta



### PLASTIC MICROLOOPS

**Description:** Premium Quality sterile plastic inoculating loops and needles. Precision moulded plastic microloops and needles are available in individual peel pouches, or in convenient tamper evident resealable packs of 20.

**Sizes available:** 1µl, 5µl and 10µl, plus needle version.

**Opening the bags:** Tear off the perforated section, initially try to gently pull the two layers apart to open the Ziplock, Note: if the package does not open easily, slide the Ziploc horizontally between two fingers to open the package.

**Method of use:** Remove loop from sterile plastic bag, and collect test material such as urine specimen, or bacterial colonies from an agar plate. Once used dispose of loop as clinical waste. Re seal the bags in the normal fashion by sliding fingers across the ziplock.

**Storage: No specific Storage Temperatures**

**Expiry Date: This product has a five year expiry date.**



Medical Wire & Equipment Co (bath) Ltd  
Corsham  
Wiltshire SN13 9RT UK  
[info@mwe.co.uk](mailto:info@mwe.co.uk) [www.mwe.co.uk](http://www.mwe.co.uk)



Revision: 4 Date 10/03/2020



Advena Ltd, Tower Business Centre, 2<sup>nd</sup> Fl.,  
Tower Street, Swatar, BKR 4013, Malta