



LyfoCults[®] Plus

All-in-one QC microorganism System



Self-Contained System

- Comprehensive selection
- Self-hydrating
- Built-in "Streaker"
- Streak for isolation
- No glass ampoule
- Available in 2-packs

Because LyfoCults[®] Plus contains a built-in inoculation device, you can inoculate either tubed or plated media and streak the dish for isolation!



LyfoCults[®] Plus are superior to other products currently on the market. Many other products contain a heavy gel content which requires added manipulation to dissolve. LyfoCults[®] Plus are lyophilized and hydrate instantly.

No need for added:

- Sterile forceps
- Pipettes
- Sterile vials
- Hydration fluid

To purchase LyfoCults[®] Plus or for further information, please contact us.

KEY SCIENTIFIC PRODUCTS

1402 D CHISHOM TRAIL

ROUND ROCK, TEXAS 79548

VOICE 800-843-1539 FAX 888-440-4208

[click to see list](#)

www.keyscientific.com/org.htm

ATCC Licensed
Derivative