

Gram-positive, Non-sporeforming Bacilli (PAGE 1 OF 2)

K1464
NAG- N-acetyl-b-d-glucoaminidase
MBGL--4mu B glucosidase
PRO-Proline naphthylamide

K1263
AGLU-Alpha glucosidase
MAGA-4mu alpha galactosidase
LGY-Leucyl-Glycine a-naphthylamide

Supplemental tests:
 K1490 WEETAB ONPG
 K2375 Spot indole

NAG

Positive
MBGL

Negative
SEE PAGE 2

Positive
PRO

Negative
AGLU

Positive

Negative

L. acidophilus
L. casei

L. catenaforme

Positive
LGY

Negative
LGY

Positive

Negative

Prop. acnes (Ind+)
Actinomyces meyeri (Ind -)

E. aerofaciens

Positive

Negative

Actinomyces bovis
Actinomyces pyogenes (ONPG+ 86%)
Prop. Acnes (Ind + 88%)

E. aerofaciens

Gram-positive, Non-sporeforming Bacilli (PAGE 2 OF 2)

K1464

NAG- N-acetyl-b-d-glucoaminidase

MBGL--4mu B glucosidase

PRO-Proline naphthlyamide

K1263

AGLU-Alpha glucosidase

MAGA-4mu alpha galactosidase

LGY-Leucyl-Glycine a-naphthlyamide(

Supplemental test

K1650 WEETAB Urease

