



The sky's the limit...

you design it ... any size ... any color ... any configuration

Pre-Printed Plates & Microscope Slides

* affordable * small minimums * minimum lead-time

Choose from one of our patterns or design your own


These specially cured standard hydrophobic or "Super Hydrophobic" ceramic inks resist harsh laboratory chemicals. These microscope slides are of the highest quality. Slides available with plain or frosted ends, with ground glass edges and 90° corners. Choose from either standard hydrophobic or "Super Hydrophobic" ink in red, green, blue, black or white. "Super Hydrophobic" inks are ideal for growing cells on the slide. Slides using hydrophobic inks can be used for standard fluorescent antibody techniques. Standard hydrophobic inks come in matte finish. FREE artwork and no screen charges. Sold by the gross (144 slides). Custom orders cannot be returned.

MINIMUM ORDER OF 5 GROSS ON ALL PRINTED MICROSCOPE SLIDES.

- 041 Custom Printed Slides (epoxy) 5-19 gross
20+ gross
- 041-T Custom Printed Slides with Super Hydrophobic Ink
- 074 Custom Printed Plates on SPECIAL SIZES INQUIRE
- 045 Plain Blank Slides (case 10 gross)
- 050 Frosted Blank Slides (case 10 gross)


Please use the following checklist when ordering:


1. Specify Catalog Number
2. Specify Quantity
3. Standard Hydrophobic or Super Hydrophobic Slides?
4. What color ink do you desire?
5. Plain, Frosted or Ceramic end?



NEW

NEW





Agglutination Slides

ideal for latex and other agglutination techniques

These reusable 2 inch x 3 inch printed glass slides are ideal for latex and other agglutination techniques. Each type of slide has sufficient space within the well to clearly see reactions. Ink is made of a Super Hydrophobic epoxy to contain a large amount of liquid within each well. Available by the gross. Inquire for custom designs and pricing.

074-1 Agglutination Slide-square well
074-2 Agglutination Slide-round well



Sterile Slides

These sterile slides are frosted 1"x 3" frosted slides that have been packaged in envelopes and autoclaved for those special uses that need sterile slides. Each envelope contains one slide. There are 50 envelopes per case. Printed or ceramic colored end slides are also available, please inquire for pricing.

050-S Sterile slides (50)




Printed TB Slides

ideal for making smears

These popular slides are ideal for making smears where mycobacterial disease is suspected. The large 20mm well accommodates a large amount of material. These black epoxy slides can be used for either acid fast or for fluorescent staining.

MINIMUM ORDER OF 5 GROSS.
041-0120 Printed TB Slides



Gram Stain Printed Slides

no more wax pencils

Avoid wax pencil lifting off the slide and leaving artifacts after staining. The black hydrophobic ink prevents mixing of smears. Ink resistant to most laboratory chemicals. Three patterns available.

MINIMUM OF 5 GROSS OF SLIDES.
030-1 Printed Gram Stain Slides
030-2 Printed Gram Stain Slides
030-3 Printed Gram Stain Slides



Grid Slides

easy to find specimens

This ingenious slide has a microscope grid on the reverse side of a frosted end slide. Ink is resistant to stains and common laboratory solvents. Each square grid fills one low-powered field (approx 100x magnification).

MINIMUM OF 5 GROSS OF SLIDES.
044 Grid Slides



Cyto-Spin Slides

specimens in the exact position you need them

These preprinted slides are designed to make locating the centrifuged smear easier. Each well is placed in the exact position where the centrifuged specimen is deposited. Available for several cyto-centrifuge systems. The ink is a black permanent hydrophobic epoxy that is resistant to most laboratory solvents. Inquire about adhesive coating on these Cyto-Spin Slides.

MINIMUM ORDER OF 5 GROSS.
056-1 Cyto-Spin Wescor (12mm well)
057-1 Cyto-Spin Shandon 1 well (8mm well)
058-1 Cyto-Spin Shandon 2 well (8mm wells)
060-1 Cyto-Spin StatSpin™ (16x24mm)



Adhesive Coated Slides

avoid losing valuable specimens

Poly-l-lysine or Silane coated slides are ideal for use with hybridization techniques. Each production run is put through extensive quality control. These adhesives have been shown to be effective but will not interfere with staining.

MINIMUM ORDER OF 5 GROSS.
067 Poly-l-lysine Slides
068 Silane Slides



Hybridization Slide

optimal adhesive qualities

These new and improved coated slides provide optimal adhesive qualities for performing hybridization and other staining techniques using tissue. Each custom batch of slides are quality controlled with brain tissue using full hybridization technique at various temperatures and with organic solvents. They meet our strict quality guidelines.

MINIMUM ORDER OF 5 GROSS.
063 Hybridization Slide
063P Hybridization Slide with custom printed design

Aqua-Hold pap pen *replaces expensive PAP-Pens*

The Aqua Hold pap pen draws a liquid repellent area around a section of tissue so that you do not waste expensive aqueous solutions. The Aqua Hold pap pen is used just like a paint pen to draw a well on a slide. This water repellent barrier creates the proper surface tension to hold aqueous solutions such as costly antibody and other valuable specimens. After application the barrier material is insoluble in alcohol and acetone but can be removed after staining with xylene or SDL's environmentally-safe Paraclear. Each pen holds 13ml of long lasting barrier material which will make hundreds of wells. Available individually.

9804 Aqua Hold pap pen (each)



Sterile Disposable Instrument Packs

pre-packaged and ready-to-use

These instruments are ideally suited for the busy laboratory where sterile instruments are needed at a moments notice. Avoid washing and exposure to infectious diseases by simply disposing of them after use. Each sterile instrument pack (cat #903) contains one sturdy stainless steel four inch forceps, one rigid plastic scalpel holder with a #11 blade and one stainless steel pointed scissor. Each order contains 12 individual sterile packs. Scissors and forceps can be re-sterilized by autoclaving.

- 903 Sterile Instrument Packs (12/package)
- 904 Sterile Forceps (25/package)
- 906 Sterile Instrument Pack with Dissection Board (10 packs)
- 9003 Sterile Scalpel #11 blade (100/package)

Slide Forceps *extra-wide grips allow for greater stability*

The extra wide tip of these handy forceps holds microscope slides securely while protecting your hands from getting stained. Autoclavable and made of durable plastic, the forceps last as long as comparable metal forceps at a fraction of the cost Packaged in lots of 10.

- 902 Slide Forceps (10 per package)
- 902-S Slide Forceps STERILE (100 per package)

Sterile Disposable Forceps *economically priced*

Sterile forceps (cat #904) are exceptionally economical and are available in individually wrapped sterile packets. Forceps are packaged 25 per package. Also available non-sterilized.

- 904 Sterile Forceps (25/package)
- 904-NS Non-Sterile Forceps (100/package)



Sterile Instrument Packs



Sterile Disposable Forceps



Slide Forceps



Fine-Tipped Forceps

Disposable Fine-Tipped Forceps *extremely affordable*

These disposable forceps facilitate work with smaller specimens. They are made of stainless steel and coated with an autoclavable plastic. Length is 5.25 inches. Packaged in units of twelve or fifty. Also available in sterile wrap, inquire for pricing.

- 9004 Fine Tipped Forceps (12/package)
- 9004-S Fine Tipped Forceps STERILE (100/package)

ThermoTote Incubator

Inexpensive variable temperature transports that plug into your car's cigarette lighter



The ThermoTote family are inexpensive, yet reliable, incubators/transporters that will operate on 12 volts either in your laboratory or alternatively in your automobile by plugging into the cigarette lighter. They take only a few minutes to come up to temperature. They also double as extra incubators for special temperature studies.

Purchase the 120V AC/DC Converter (Cat. No.7155) to allow ThermoTote to easily plug into a normal power outlet. Call to inquire about ThermoTotes with self-contained batteries. Also call to inquire about customized ThermoTotes or Service Contracts.

BATTERY POWER (for use with Small ThermoTote12) When power is not available. Can be recharged from 12 volt or 120V source (optional converter). Trickle mode allows the battery, once charged, to deliver the maximum power necessary. Included is 14.4 volt battery pack and charger. Call to inquire about Batteries for other ThermoTote sizes.



TEMP RANGE	Small 25°C to 50°C HEAT only	Medium 25°C to 50°C HEAT only	Deluxe 2°C to 65°C HEAT or COOL
POWER	12 volts	12 volts	12 volts
INTERNAL DIMENSIONS	7.75"x10.75"x7"	7.75x10.75"x13.25"	17"x12"x13"
EXTERNAL DIMENSIONS	13" x 9.5" x 8.75"	13" x 9.5" x 15"	22" x 15.25" x 16.75"
WEIGHT	4 pounds	4.75 pounds	17 pounds
INCLUDES	cigarette lighter power adapter plug	cigarette lighter power adapter plug	cigarette lighter power adapter plug & 1 shelf included

- 7145 Small-Size ThermoTote12 Adjustable with LCD Display
- 7148 Small-Size ThermoTote12 Adjustable with LCD Display & BATTERY
- 7146 Medium-Size ThermoTote24 Adjustable with LCD Display
- 7151 ThermoTote Deluxe Adjustable

- ACCESSORIES for use with ALL ThermoTote products unless noted
- 7155INTL International Converter 220v/240v
 - 7155USA Plug 120v (for standard U.S. plugs)
Converts 120 volt power to 12 volt DC for non-portable use
 - 7160 Cigarette Lighter Adapter (use with only Large "Deluxe Incubator")
 - 7162 Cigarette Lighter Adapter (use with "Small" or "Medium Incubator")
 - 7157-12 Wire Petri Dish Racks 4 racks
holds up to 10 plates per rack (use with only "Small Incubator")
 - 7157-24 Wire Petri Dish Racks 4 racks
holds up to 10 plates per rack (use with only "Medium Incubator")
 - 7152-01 Multiple Plug allows for three units to function at the same time



Perfect for courier drop off and pick-up



Plugs into car cigarette lighter

Snap n' Digest Sputum Digestion System

*Pre-measured, ready-to-use
less expensive than the competition*

Each bottle contains 75ml or 150ml of an alkali solution with a sealed glass ampoule of preweighed N-acetyl-L-cysteine (NALC). When needed, the glass ampoule snaps easily with two fingers and then the contents rapidly dissolve in the fluid. The unbroken capsule is stable for two years without refrigeration. Each bottle has a specially designed tube to help you easily trap and break the NALC capsule. With no weighing of components and no wasted chemicals, this self-contained system will save you time and money. Included in the kit is sufficient dry buffer powder needed to complete the neutralization of the digestion process. This inexpensive digestion fluid for processing mycobacterial specimens is easy to use and avoids time spent weighing reagents. Each package contains 10 non-breakable bottles and 5 packets of dry TB buffer powder for washing. SDL also offers a larger size Snap n' Digest™ of 150ml fill. For those laboratories with higher volume, this larger size will save you both time and money.



TB Buffer, Pre-Sterilized

ready-to-use and economical

Pre-Sterilized Buffer consists of economical bottles containing ready-to-use TB buffer for washing sputums after digestion and decontamination. The buffer has a long expiration date. Available in 100ml or 500ml bottles. Also available in non-sterilized powdered form.

669-100 Sterilized TB Buffer (10 x 100ml)

669-500 Sterilized TB Buffer (500ml)

6671 Powdered Buffer Packets
(5 x 4.2g each)

667 Snap N' Digest™ (75ml)

668 Snap N' Digest™ (150ml)



Slide Heater *safely fix, dry & stain slides*



This low-profile slide heater is uniquely safe because it runs on a low 12 volts. Measuring 10" x 4" it is ideal to fix, dry, or stain up to 20 slides. Thermostatically controlled, the temperature can be set anywhere from room temperature to 65°C. This large heater holds up to twenty (20) slides. Hot surface is made of a resin powder coated aluminum.

504-20 Slide Heater

Mycobacteria Inactivation Fluid & Rack

In a recent paper¹, a rapid method was described to protect workers from aerosolization of tubercula bacilli on slides prior to removal from the bio-safety cabinet. Treatment with a disinfectant aided in fixing and inactivating potentially pathogenic organisms. This reagent is now available in a 500ml glass bottle with an accompanying, convenient and easy-to-handle, 125ml spout bottle for each application. To facilitate disinfection in biohazard hoods, a low-profile staining rack and dish is also available. Polyethylene Staining Dish and Rack is resistant to activity of phenol. Rack holds up to seven slides.

188 Mycobacteria Inactivating Solution (500ml & 1 squirt bottle)

189 Mycobacteria Inactivation Staining Dish with Rack

¹ J.Clin.Micro.42:3, 1288.

AcriFluor™

rapid detection method for Mycobacteria, as well as other acid fast organisms

This new fluorescent technique for documentation of Acid Fast organisms in tissue and smears is both rapid and specific. The new modified AcriFluor stains all acid fast organisms so that you can find them easily. For Mycobacteria in sputum, the organisms simply light up on the slide. This non-antibody technique is extremely cost effective. The procedure takes less than twenty minutes to perform. One kit will enable you to do hundreds of assays. The kit comes with one 125 ml bottle of AcriFluor™ and one 125 ml bottle of Decolorizer.

384 AcriFluor Stain Kit (125ml)

TB Splash-Proof Containers


maximum safety in discarding contaminated liquid waste

These containers are made according to the CDC's recommendations to give maximum safety in discarding contaminated liquid waste. Entire polypropylene container and funnel assembly is autoclavable and reusable. Each set includes two containers and one tip-proof acrylic rack.

251 TB Splash-Proof Container (500ml capacity)

252 TB Splash-Proof Container (1000ml capacity)

252-999 Funnel & Stopper Set (4 replacements)



Acid Fast Control Slides, Kinyoun Stains, & Ziehl Neelson Stain

convenient control slides & stains

These unique Acid Fast Control Slides simulate a clinical sputum smear with appropriate background material. Slides have two wells containing a positive and negative smear. Can be used with Kinyoun (Cat. No. 353), Ziehl-Neelsen or SDL's new product AcriFluor™.

- 350 Acid Fast Stain Control Slides (10)
- 351 Acid Fast Stain Control Slides (50)
- 353 AFB Stain/Kinyoun Stain Set (250ml)
- 353-A AFB Kinyoun Stain (3x250ml)
- 353-B AFB Kinyoun Decolorizer (3x250ml)
- 353-C AFB Kinyoun Counterstain (3x250ml)
- 358 Acid Fast Tissue Control Slides (10)
- 377 Ziehl-Neelson Modified AFB Kit
(Stain & Brilliant Green, 250ml each)

Printed TB Slides

ideal for making smears

These popular TB slides are ideal for making smears where mycobacterial disease is suspected. The large 20mm well accommodates a large amount of material. These black epoxy slides can be used for either acid fast or for fluorescent staining. MINIMUM ORDER OF 5 GROSS.

041-0120 TB Slides



Trehalose (RAT) Test

rapidly distinguishing *C. glabrata*

Studies¹ show that SDL's new RAT assay can rapidly and affordably distinguish *Candida glabrata* from other yeasts. Modeled after Mayo Clinic's procedure, this new rapid assay system allows for easy identification of *C. glabrata*. Each slide contains four wells of substrate, you can use them all at the same time or each well individually. Simply add one drop of reconstituting fluid (included in kit) and mix in organism. Read color reaction within two (2) hours. Each package contains enough for 40 assays.

¹ T. Scognamiglio, L.P. Kiss, D.H. Larone, ASM poster C-105:164, 2004.

761 Trehalose Slides (10 slides)

Evaluation of a New Commercially Available Rapid Assimilation of Trehalose (RAT) Test for the Identification of *Candida glabrata*

T. Scognamiglio, L.P. Kiss, D. H. Larone,
Cornell Medical Center, NewYork-Presbyterian Hospital, New York, NY



Conclusion: The new SDL RAT assay is very useful for distinguishing *C. glabrata* from other *Candida* species. Optimal results are obtained with isolates from BD's CAC (sensitivity & specificity 100%). One hour incubation is best, as isolates from SDA may give false positives if incubated for 4 hours. Additionally, isolates of *C. glabrata* from SDA may give false negative results and require further testing for identification. Of note, the SDL assay closely mimics the previously reported results of the NCCLS (M35-P) proposed method and has a more rapid reaction time and is more reasonably priced than other products currently on the market.

ASM Poster 2004, C-105



India Ink Capsule Control

Now have a standardized permanent, preservative fixed control for use with India Ink. This suspension of *Cryptococcus* is grown to maximize capsule production. Each vial contains 1ml of fluid containing many organisms. India Ink Negative Control is a 1ml suspension of *Candida albicans*.

- 756 India Ink Capsule Positive Control(1 ml)
- 755 India Ink Capsule Positive Control (1ml) with 6 x 3.5ml india ink
- 756-1 India Ink Negative Control(1ml)
- 754 India Ink (6 x 3.5ml bottles)



Aluminum Mycology Needles

Unlike wooden-handle needles, these sturdy aluminum mycology needles will not burn or melt. They come in a set of one straight and one curved.

- 687 Aluminum Mycology Needles (one each curved & straight)

NEW

Methylene Blue/ Glucose Enhancement Reagent

Methylene Blue/Glucose reagent (MBG) is currently used to aide in reading clearer zones of fungal inhibition when performing disc susceptibility tests. They are available in two forms. The contents of our Kwik Pour MBG tubes can be melted in boiling water and poured into Petri Dishes as needed or the Methylene Blue/Glucose reagent can be flooded onto pre-made Mueller Hinton Plates. Kwik Pour Tubes have a 1 year shelf life. Each tubes contains enough media for one plate. There are 25 tubes per box. The 1.5 ml of MBG reagent is sealed in easy to snip pipettes. Twelve pipettes per box.

- 650 Kwik Pour MBG Tubes (25 tubes/box)
- 651 MGB Reagent (25 pipettes/box)

NEW



KOH-Methylene Fungal Control Suspension

Now you can quality control the Methylene blue-KOH. This 1ml suspension is filled with inactivated mycelia and spores from an *Aspergillus* sps. along with other cellular debris. Makes an ideal control for a skin scraping.

- 6025 KOH Fungal Control (1ml vial)

Mycology Teaching Slides

beautiful, quality images

Make your own sets of mycology teaching slides with these high-quality glass microscope slides. All fungi are stained and mounted with SDL's MycoPerm. Minimum of any ten slides per package.



- 670 thru 675 Mycology Teaching Slides (Minimum of any ten slides per package, each slide)

- | | | |
|-------------------------------------------|----------------------------------------|-------------------------------------|
| 670-02 <i>Streptomyces griseus</i> | 673-05 <i>Curvularia</i> sp. | 674-08 <i>M. vanbreuseghemii</i> |
| 670-03 <i>Nocardia</i> sp. | 673-06 <i>Drechslera</i> sp. | 674-09 <i>T. ajelloi</i> |
| 671-01 <i>Blastoschizomyces capitatus</i> | 673-10 <i>Nigrospora</i> sp. | 674-11 <i>T. equinum</i> |
| 671-02 <i>Candida albicans</i> | 673-12 <i>Phialophora verrucosa</i> | 674-12 <i>T. megninii</i> |
| 671-03 <i>Cryptococcus neoformans</i> | 673-15 <i>Scedosporium apiospermum</i> | 674-13 <i>T. mentagrophytes</i> |
| 671-04 <i>Geotrichum candidum</i> | 673-16 <i>Scedosporium prolificans</i> | 674-14 <i>T. rubrum</i> |
| 671-05 <i>Hansenula</i> sp. | 673-18 <i>Ulocladial</i> sp. | 674-15 <i>T. schoenleinii</i> |
| 671-07 <i>Sporothrix schenckii</i> | 673-19 <i>Wangiella dermatitidis</i> | 674-16 <i>T. terrestris</i> |
| 672-04 <i>Mucor</i> sp. | 674-01 <i>Epidermophyton floccosum</i> | 674-17 <i>T. tonsurans</i> |
| 672-06 <i>Rhizopus</i> sp. | 674-02 <i>Microsporium audouinii</i> | 675-01 <i>Acremonium</i> sp. |
| 672-07 <i>Rhizomucor</i> sp. | 674-03 <i>M. canis</i> | 675-06 <i>Aspergillus flavus</i> |
| 672-09 <i>Syncephalastrum</i> sp. | 674-07 <i>M. gypseum</i> | 675-07 <i>Aspergillus fumigatus</i> |
| 673-01 <i>Alternaria</i> sp. | | 675-09 <i>Aspergillus niger</i> |
| 673-02 <i>Bipolaris</i> sp. | | 675-16 <i>Chrysosporium</i> sp. |
| 673-03 <i>Chaetomium</i> sp. | | 675-18 <i>Fusarium</i> sp. |
| 673-04 <i>Cladosporium</i> sp. | | 675-19 <i>Paecilomyces</i> sp. |
| | | 675-20 <i>Penicillium</i> sp. |
| | | 675-32 <i>Scopulariopsis</i> sp. |
| | | 675-33 <i>Sepedonium</i> sp. |



MycoEase

faster, safer fungal identification



This revolutionary, patented and cost effective slide system helps in the identification of fungi faster and with increased safety

Each of the two sealed slide chambers contains a cellulose cylinder with dehydrated media. To initiate, rehydrate the media with enclosed reconstituting fluid. Then inoculate from either a direct source or by subculture. One chamber contains Sabouraud Dextrose medium and the other Potato Flake medium.

Safer and works faster than traditional methods. The humidity chamber provides an ideal environment for fungal growth. The patented MycoEase allows for air exchange to optimize growth but spores cannot escape the chamber. The unit is constructed so that viewing growth of molds can be performed safely each day. Specially treated

coverslips attract spores so that your focus is optimized. It takes just seconds to put the whole unit on the microscope stage and view daily growth. No more guessing when your culture is ready. Once the identifying fungus is viewed, chamber can be discarded or permanent slide preparations can be made using MycoPerm.

SHELF LIFE IS ONE YEAR. MycoEase is ideal for field work or teaching mycology. Specially prepared coverslips, rehydrating pipettes, and materials for 30 tests included with each case.

6000 MycoEase (15 slides includes material for 30 tests)

753 MycoPerm - BLUE (12ml)

752 MycoPerm - RED (12ml)



MycoPerm BLUE & MycoPerm RED

no more nail polish, a unique permanent Mycology Stain & Mounting Fluid with acid fuchsin or lacto-phenol cotton blue

Studies(1) now show that laboratories can keep better records by retaining mycology preparations indefinitely. MycoPerm makes using messy nail polish obsolete. Tease mount preparations can be made and stored for long periods of time. The novel MycoPerm stain contains plasticizers which help in making a durable long-life preparation. Available in dropper bottles of 12ml and 100ml refill size. The **Mycology Starter Kit** gives the mycologist all of the necessary reagents to aid in the identification of fungi. Each kit comes with six (3.5ml) dropping bottles of 10% KOH, Brite-Fluor, one bottle 12ml of MycoPerm and a durable rack.

(1) Perry et al., ASM poster 125:70, 1996.

753	MycoPerm - BLUE	12ml bottle
752	MycoPerm - RED	12ml bottle
750	Mycology Starter Set BLUE	MycoPerm-BLUE, 6x3.5ml KOH, 3.5ml bottle Brite-Fluor, rack
750-9	Mycology Starter Set RED	MycoPerm-RED, 6x3.5ml KOH, 3.5ml bottle Brite-Fluor, rack
751	10% KOH	6 x 3.5ml dropping bottles
759	Bright-Fluor	6 x 3.5ml bottle

Fungi-Tape

Fungi-Tape is an efficient way of performing slide mounts of fungi from primary cultures. The adhesive tape is pliable enough for you to make a tape loop and then sample your culture. This product is fully compatible with MycoPerm to create a long lasting coverslip. Each roll contains approximately 300 coverslips.


745 Fungi-Tape (1 roll)



Mycology Images CD-Rom *fabulous for presentations or teaching*

Beautiful images of all the major fungi that can be used for web-sites, teaching or any other type of presentations. Simply insert the images into any Powerpoint presentation. All images in .jpg format. Not intended for printing purposes. Images listed below.

6790 Mycology CD-Rom



Alternaria sp.	Candida albicans	Microsporium vanbreuseghemii	Penicillium thomii	Trichophyton megninii
Aspergillus caespitosus	Chaetomium sp.	Mucor	Penicillium variable	Trichophyton mentagrophytes
Aspergillus carbonarius	Chrysosporium sp.	Nigrospora sp.	Phialophora verrucosa	Trichophyton rubrum
Aspergillus clavatus	Cladosporium sp.	Nocardia sp.	Rhizomucor sp.	Trichophyton schoenleinii
Aspergillus flavus	Coccidioides immitis	Paecilomyces sp.	Rhizopus sp.	Trichophyton tonsurans
Aspergillus fumigatus	Cryptococcus	Penicillium citrinum	Scedosporium	Wangiella
Aspergillus japonicus	Curvularia sp.	Penicillium chrysogenum	Scedosporium prolificans	
Aspergillus niger	Drechslera sp.	Penicillium chrysogenum	Scopulariopsis sp.	
Aspergillus parasiticus	Epidermophyton floccosum	Penicillium expansum	Sepedonium sp.	
Aspergillus sojae	Fusarium sp.	Penicillium funiculosum	Sporothrix schenckii	
Aspergillus tamarii	Geotrichum	Penicillium chrysogenum	Streptomyces	
Aspergillus terreus	Hansenula sp.	Penicillium janczewskii	Syncephalastrum sp.	
Aspergillus versicolor	Microsporium audouinii	Penicillium restrictum	Trichophyton ajelloi	
Bipolaris sp.	Microsporium canis	Penicillium roquefortii	Trichophyton equinum	
Blastoschizomyces capitatus	Microsporium gypseum	Penicillium simplicissimum		

Scientific Device Laboratory

Order from: KEY SCIENTIFIC PRODUCTS 800-843-1539

MycoEASE 7 page

Shrink Seals

safe & secure sealing of Petri Dishes




maintain a safe, tight seal and keep agar plates longer by retaining moisture

Shrink Seals -- are colorless transparent cellulose bands for sealing Petri dishes and test tubes. They make an invaluable addition to any microbiology laboratory. Fantastic for in-house prepared plates that are not used quickly. Ideal for fungus cultures that require prolonged incubation and for sensitivity plates for TB. Easy to apply. The size of the standard band is 92mm x 25mm. Test tube sizes listed below. Now available in containers of 400 seals.

Stacking Blocks -- are available for those labs with limited table top space. This spacing device allows proper space between your plates while they are drying and sealing.

Test Tube Shrink Seals -- Keeps a secure seal on test tube caps. Easy to apply. Test tube sizes listed.



261 Shrink Seals 1 container, 400 seals

259	Stacking Blocks	18 blocks
266	Shrink Seals Test Tubes- 60mm x 25mm	3 gross
267	Shrink Seal Test Tubes- 16mm x 15mm	3 gross
268	Shrink Seal Test Tubes - 18mm x 15mm	3 gross
269	Shrink Seal Test Tubes- 21mm x 20mm	3 gross
270	Shrink Seal Test Tubes- 27mm x 25mm	3 gross
271	Shrink Seal Test Tubes- 35mm x 25mm	3 gross



Cryocare Bead Storage System

long-term storage of bacteria, fungi, yeasts, QC cultures, cultures, or any other organisms

The bead system can even store fastidious or delicate organisms. Inoculate the beads, shake, draw off liquid with sterile pipette, then freeze. Viability at -70°C can be up to 10 years, 1 to 3 years at -10°C to -15°C. Simply remove a bead and run on agar plate or drop into broth media to grow. Beads packaged 20-25 beads per vial, 50 vials per box. Color code organisms by purchasing various colors.

110 CryocareBead (50 vials/box) Specify color when ordering by adding suffix: red (-01), blue (-02), yellow (-03), green (-04), white (-05), or mixed colors (-06)

Super-Chlor™ Bleach Pad prediluted bleach wipes

These prediluted bleach wipes are ideal for decontaminating environmental surfaces and equipment. Protects against both viruses and bacteria. Each towelette is impregnated with sodium hypochlorite, precisely prediluted with distilled water and sealed in a gas impermeable easy-open packet. The standard size for each SuperChlor™ is 9" x 9".

9817 SuperChlor™(100 packets)

Stainless Steel Anaerobic Jar

affordable incubators

These unique and inexpensive anaerobic jars are perfect for both the small and large laboratory. These jars seal easily and hold either an anaerobic environment or carbon dioxide. Our sturdy rack fits comfortably into our anaerobic jars and allows the user to withdraw all of the plates at one time. The small jar holds eight plates with our rack or nine plates without. The large jar holds eleven plates with the rack, twelve without. All polished stainless steel jars have clear acrylic lid with clip to view Methylene blue test strip indicator to measure anaerobiosis. Sufficient room is left for the anaerobic envelope and Methylene blue detector system.

714	Small Stainless Steel Jar - 9 plates no rack included
715	Large Stainless Steel Jar - 12 plates no rack included
716	Small Stainless Steel Jar with rack - 8 plates
717	Large Stainless Steel Jar with rack - 11 plates
722-S	Replacement Seals for Stainless Steel Incubator (5)

Anaerogen™ Fast Anaerobic Atmosphere Generation System

AnaeroGen™ enables oxygen to be safely consumed with no hazardous substances, no hydrogen by-product and no dangerous buildup of pressure in the jar without the help of catalysts or water. The final atmosphere is less than 1% oxygen, supplemented with carbon dioxide - ideal for optimum growth of anaerobes. Ten envelopes per box. Indicator Test Strips verify that the atmosphere, once activated is truly anaerobic. Packaged in a box of 100 strips.


725	Anaerogen™ Pouches (10/box)
726	Anaerogen™ Indicator Strips (100)



QC Boards the neatest way to keep QC records

These Quality Control Boards are specially designed to record and view laboratory data. Keep your CAP or JCAH quality control records clean and intact. Take board off surface, write your QC data and replace, with no cumbersome handling of books. Boards cling to either iron or non-ferrous surfaces (walls, doors or sides of machinery). Data is protected by a clear mylar cover. Backing is made of 1/8" acrylic which measures 8.5" x 11". Packaged in lots of six. Not recommended for doors that are forcibly closed.

100	QC Board Steel Surface (6 boards/pkg)
101	QC Board Non-Steel Surface (6 boards/pkg)



Teaching & Reference-Atlas of Parasitology CD-ROM

This unique, user-friendly CD-ROM of the *Color Atlas of Parasitology* is an ideal teaching tool and reference for all laboratories to have in their library. It contains a comprehensive collection of photographs of sections, life cycles, wet mounts and case histories for most of the human parasites seen in the clinical laboratory. Each of the photos and diagrams have arrows for important areas to view. Windows™ 95, 98, NT & Macintosh™ System * or above, requires Web browser (HTML 2.1 or above compliant).



025 Atlas of Parasitology CD-ROM



Parasafe™ formalin-free fixative & Mini-Parasafe™

Safe... Parasafe® is nontoxic and non-carcinogenic. This rapid fixative meets OSHA requirements for safety and disposal. Parasafe® contains no formaldehyde or mercury.

Fast... The patented Parasafe® procedure cuts stool processing time to under one hour. Slide preparation dries rapidly. Improved rapid Trichrome stain procedure saves time. Eliminates the iodine-alcohol step for staining thus saving additional time. Our patented Slide Fixative assists in fixing parasites to slide for staining.

Economical... Parasafe® is half the cost of the two vial system.

Effective... Parasafe® allows for better patient compliance. Now you can have a stool fixative for wet mounts and for permanent stains in one tube. The concentrated stool specimen can be analyzed by wet mount and Trichrome stain at half the price of formalin and PVA. Results give optimum definition of ova and fine protozoal nuclear detail. No shrinkage.

Other Options--Mini-Parasafe is Parasafe® packaged in 15ml vials. Both standard and Mini-Parasafe are compatible with our Fecal Parasite Concentrators (Cat# 6052). Also available is Cary-Blair, a medium for culture of feces. SAF is an inexpensive fixative used for concentration of feces and making permanent smears using either Iron Hematoxylin or Trichrome Stains. Available in packages of 24 tubes with a 15ml fill or "Mini" size with 100 tubes with a 7.5ml fill.


Parasafe® is a patented formalin-free fecal fixative which combines all the qualities of formalin and PVA into one tube without toxic effects

Parasafe® is HALF THE COST of the other two vial systems.

Unlike formalin, Parasafe® is a nontoxic and non-carcinogenic stool fixative for wet mounts and for permanent stains in one tube

Ordering & Pricing Information -- In addition to filled Parasafe® vials and Mini-Parasafe vials, Clean Vials are available for transport of stool specimens meant for culture. Both Parasafe® vials and Mini-Parasafe are packaged in vials that are compatible with our closed system concentrators. Parasafe® & Clean Vial Set with Sealable Transport Bag includes: 1 vial of Parasafe®, 1 Clean Vial and a Sealable Transport Bag (with protective pouch for specimen paperwork) in order to make specimen collection safer, easier and more convenient. Packaged 12 sets per order.

6030	Parasafe® (24 vials)
6031	Parasafe® & Clean Vial Set (12)
6038	Clean Vials (24)
6046-100	Mini-Parasafe (100)
6032	Slide Fixative (7.5ml)
6036	Carey Blair only(24)
6065-100	SAF Mini Vials (100)
6060	SAF Vials 24 vials



Paraclear® Environmentally-Safe Clearing Agent

PARACLEAR® avoids all the problems of using and disposing of xylene. Paraclear® is a xylene-free clearing agent. Unlike xylene or citrus limonene, it has no overpowering odor. Use Paraclear® without any of the health concerns of xylene. This xylene replacement gives you excellent clinical results. Fast evaporating Paraclear® is ideal for use with all Trichrome procedures. It gives consistent results with excellent clarity. Paraclear® effectively deparaffinizes and eliminates residual chemicals. Paraclear is convenient and is less expensive to use. Paraclear® is completely soluble in alcohol and has no effect on most laboratory plastics. It is also less flammable than xylene. Costly distillation becomes unnecessary.

None of the disposal headaches of xylene... none of the odor of D-limonene

NO EXTRA HAZARD SHIPPING COSTS are associated with this product.

- 615-1 Paraclear® (1 gallon)
- 615-4 Paraclear® (4x1 gallon)
- 615-5 Paraclear® (5 gallon cube)



Fecal Concentrator Kits *cost effective kits*

This cost-effective disposable Fecal Concentrator Kit contains all the items needed for concentration of feces for ova and parasites. The kit contains 50 each of standardized fine mesh filter funnels, 15ml polypropylene centrifuge tubes, tube caps, one 15ml bottle of Mucolytic Digesting Fluid and Parawash (250ml). Kits are available without Parawash for use with our Parasafe fixative. Items also sold separately.


Closed concentration system (cat # 6052 & 6056) allows for increased safety with a minimum of handling and odor. The 0.6mm holes in the filter allows the parasites to pass through but blocks larger fecal material. Available in 15ml tubes or 50ml Jumbo (holds 30ml) collection tubes. Each package includes a preassembled tube-strainer unit, Fecal Mucolytic reagent, dispensing spoons (for 15ml tubes), processing trays, and applicator sticks. Pictured with SDL's environmentally-safe, formalin-free fixative Parasafe®.

- 641 Fecal Concentrator Kit (50 funnels/tubes)
- 642 Fine Mesh Filter Funnels (250)
- 643 Centrifuge Tubes & Caps (250)
- 644 Mucolytic Digesting Agent (15ml)
- 052 Closed Fecal Parasite Concentrators (15ml,100/box)
- 6056 Closed Fecal Parasite Concentrators (50ml,100/box)



Order from: KEY SCIENTIFIC PRODUCTS 800-843-1539

Acrifluor™ *rapid detection fluorescent stain for*



This new fluorescent technique for documentation of Acid Fast organisms in tissue and smears is both rapid and specific. The new modified AcriFluor™ stains all acid fast organisms so that you can find them easily. Whether you have Cryptosporidia, Nocardia & Isospora in stool, Nocardia in tissue, or Mycobacteria in sputum, the organisms simply light up on the slide. A single organism is easy to find. This non-antibody technique is extremely cost effective. The procedure takes less than twenty minutes to perform. One kit will enable you to do hundreds of assays. The kit comes with one 125 ml bottle of AcriFluor™ and one 125 ml bottle of Decolorizer. SDL's printed TB slides make a valuable addition to the stain.

384 AcriFluor™ Stain Kit (125ml)
041 -0120 TB Slides

DeToX *Environmentally-Safe Formalin Neutralizer*

Dispose of formalin without safety concerns, special storage or handling charges. Once specimen is ready for disposal, pour the DeToX into the container, shake, and within minutes the formalin is fully neutralized and can be disposed of safely. DeToX is an EPA approved waste modification reagent. One test reagent is included with each order to assure complete neutralization.

617 DeToX Neutralizer (one 10lb. Jar)
617-100 DeToX Neutralizer (100 small vials)
617-99 DeToX Confirmatory Test



Zinc Sulphate Solution

ready to use solution in a variety of specific gravities

Avoid time consuming preparation of this concentration solution for ova and parasite procedures. Solutions are carefully prepared for optimal results. Zinc sulfate solutions available for either fresh preparations with a specific gravity of 1.18 or formalin preserved stool with a specific gravity of 1.20. Packaged in bottles of 250ml.

731 Zinc Sulfate Solution spec grav 1.18 (250ml)
732 Zinc Sulfate Solution spec grav 1.18 (6 x 250ml)
733 Zinc Sulfate Solution spec grav 1.20 (250ml)
734 Zinc Sulfate Solution spec grav 1.20 (6 x 250ml)



Snap n' Stain for Parasite Wet Mounts

When you don't have time to make reagents, here is a fast and efficient way to stain parasites in wet mounts

Ideal for the busy laboratory that does not have the time to make these unstable reagents. Once the ampoule in each dropping bottle is crushed, the staining solution is ready to use. The iodine preparation is made according to the formula of Dobell and O'Connor. Shelf life prior to crushing ampoule is eighteen months. Each dropping bottle contains a total of 6ml of iodine solution. Packaging contains six dropping bottles.

785 Snap n' Stain Iodine (6 bottles)



Trichrome Stain *conveniently packaged*

These stains are affordable and convenient. Packaged in 250ml easy to handle bottles.

782 Trichrome Stain(250ml)
783 Trichrome Stain (6 x 250ml)



Pinworm Collectors

paddle can be viewed directly on microscope stage

New and improved design and new adhesive make this Pinworm Collector more flexible. The clearer and thinner plastic paddle makes it easier to see ova. These Pinworm Collectors contain a paddle with one "sticky side" for pinworm collection. Less artifacts are observed. The paddle can be viewed directly on the microscope stage by snapping off the top. Instructions included. Packaged 25 vials per box.

582 Pinworm Collector (25 vials)

Explosion Proof Stool Containers

the most convenient & safe collection container

These easy-to-use stool containers will allow gases from stool specimens to escape slowly, avoiding explosion. The screw on top lets you open the container slowly and safely. The leakproof container has a wide 3.5" mouth. The container is made of non-cracking, white polypropylene. Packaged in 100 containers per box.

581 Stool Containers (100)

Gram Stain Sets *convenient & affordable*

This Gram Stain Set contains all the ingredients needed to perform a rapid and accurate stain. Bottles are easy to grip and come with pop-up spout. Set contains one 250ml bottle each of Gentian Violet, Iodine, Decolorizer and Safranin (#392-250).

392-250 Gram Stain Set of 4 bottles
392-60 Gentian Violet Six Pack (6 x 250ml each)
392-70 Gram Iodine Six Pack (6 x 250ml each)
392-80 Gram Decolorizer Six Pack (6 x 250ml each)
392-90 Gram Safranin Six Pack (6 x 250ml each)

Parawash® *a Non-Toxic alternative for ethyl acetate & Hemo-De*

Parawash® has been shown to be an effective and safe substitute wash fluid for ethyl acetate in the fecal analysis of ova and parasites when formalin is used as a preservative. Results¹ show that nontoxic Parawash gives two to three times higher recovery rate than other washes. It has no odor and is less flammable than other reagents. NOT for use with Parasafe®.

¹ Matthews et. al., ASM Ses 43:48, 1996.

6121 Parawash® (6 bottles, 250ml)

Parasitology Teaching CD-Rom

Beautiful images of all the major parasites that can be used for websites, teaching or any other type of presentations. Simply insert the images into any Powerpoint presentation. All images in .jpg format. Not intended for printing purposes. Images listed.

- 7134 Set A (22 images)
- 7135 Set B (19 images)
- 7136 Set C (17 images)
- 7137 Set A+B+C (58 images)

SET A

Entamoeba histolytica cyst
Entamoeba histolytica troph
E. hartmanni cyst
E. hartmanni troph
E. coli cyst
E. coli troph
E. nana cyst
E. nana troph
Iodamoeba butschlii cyst
Iodamoeba butschlii troph
Giardia lamblia cyst
Giardia lamblia troph
Dientamoeba troph
Chilomastix cyst
Chilomastix troph
Trichomonas troph
Blastocystis
Acanthamoeba
Cryptosporidium sps.
Cyclospora sps.
Isospora sps.
Microsporidium sps.

SET B

Plasmodium vivax ring
Plasmodium vivax troph
Plasmodium vivax schizont
Plasmodium vivax gametocyte
Plasmodium malariae ring
Plasmodium malariae troph
Plasmodium malariae schizont
Plasmodium malariae gametocyte
Plasmodium falciparum ring
Plasmodium falciparum troph

SET C (OVA)

Enterobius vermicularis
Ascaris lumbricoides
Trichuris trichiura
Necator Americanus
Strongyloides larvae
Diphyllobothrium latum
Taenia sps.
Hymenolepis nana
Hymenolepis diminuta
Dipylidium caninum
Echinococcus(hydatid sand)
Fasciola hepatica
Clonorchis sinensis
Paragonimus westermani
Schistosoma mansoni
Schistosoma japonicum
Schistosoma haematobium

NEW



Low Power Stain Heater or Water Bath

This unique 12 volt high temperature bath considerably speeds up staining of larger amount of slides. It can be used for either Trichrome or Modified Trichrome procedures. It is made of easy to clean 20 gauge stainless steel. Reaches temperatures of up to 70°C with an easy-to-read and set LCD display. Comes with rack that holds up to 36 slides. Microwave resistant racks can be purchased separately. The 120V transformer that powers the unit is UL and SA approved.

- 507 Stain or Water Bath
- 9803 Microwave Rack



NEW

StainQuick™ Trichrome

complete staining in less than 2 minutes



- 795 StainQuick™ for Trichrome Staining Kit
- 796 Replacement StainQuick™ Trichrome Stain and Decolorizer Set (250ml each)
- 792 Iodine Replacements (20 assemblies)

This rapid staining procedure gives excellent results and takes as little as 2 minutes to perform with PVA smears and less with other fixatives. StainQuick™ maintains a constant temperature of the stain. The heated stain penetrates organisms in less than a minute. Characteristic morphology of the protozoa is readily seen. The six staining 60ml coplin jars and durable rack are made of microwavable polyethylene. The system runs on 7.5 volts and components are SA and UL approved. Kit includes StainQuick™ Trichrome Stain, rack, jars, heater, iodine snap-its & decolorizer.

StainQuick™ Microsporidia

complete staining in less than 1-minute



Stain microsporidia sps. faster and more efficiently.¹ StainQuick™ maintains a constant temperature in the stain. The heated stain penetrates organisms more efficiently and in less than a minute.² The decolorizer/counterstain stains all background organisms for greater specificity. Durable 60ml coplin jars and rack are made of microwavable polyethylene. The system runs on 7.5 volts and components are SA and UL approved. For staining of microsporidia use our modified Trichrome Stain Kit² (Cat # 7843) which includes one bottle each of stain and decolorizer/counterstain (250ml each).

¹ Patten, et. al., ASM, Ses 43:40, 1996
² Patten and Lipton, ASM, Ses 94: C220, 1997


- 1951 Microsporidia StainQuick™ System
- 7843 Modified Trichrome Stain and Counterstain Set (250ml each)
- 7842M Decolorizer/Counterstain Replacement (250ml)

FlorMount

aqueous mounting media for water based reagents

This mounting medium is used when slides are processed using aqueous reagents. This stain produces no fluorescent background. FlorMount dries quickly.

- 618-100 FlorMount (100ml bottle)
- 618-400 FlorMount (4x100ml bottle)
- 618-500 FlorMount (500ml bottle)



Suremount Safe Mounting Media

xylene-free mounting medium dries clear

This new rapid drying and economical, environmentally safe xylene-free mounting medium dries clear to optimize your tissues and smears. Easy to apply. Comes in 100ml and 500ml bottles.

- 616-100 Sure-Mount (100ml)
- 616-400 Sure-Mount (4x100ml)
- 616-500 Sure-Mount (500ml)

Order from: KEY SCIENTIFIC PRODUCTS 800-843-1539

Parasitology 11 page

Stained Tissue and Blood Parasites

Each tissue smear or section contains numerous organisms and are ideal teaching aids. All slides are pre-stained and ready to be viewed. Each package contains five (5) slides of five (5) slides.

Protozoa

- 491 Plasmodium vivax (Giemsa)
- 492 Plasmodium malariae (Giemsa)
- 493 Plasmodium falciparum (Giemsa)
- 4936 Babesia sps. (Giemsa)
- 494 Trypanosoma rhodesiense
- 8967 Toxoplasma gondii (brain, H&E)

Helminths


- 553 Schistosoma mansoni (liver, H&E)
- 550 Trichinella spiralis (muscle, H&E)
- 469 Dirofilaria immitis microfilaria (bld., Giemsa)
- 496 Loa loa
- 497 Brugia malayi
- 4939 Anaplasma sp.
- 560 E. histolytica (cecum)

Bacteria

- 8925 Bacteria (Brenn & Brown)
- 8905 M. tuberculosis (Acid Fast)

Fungi (Methenamine Silver)

- 8960 Candida albicans



Protozoa Stained Slides



Every laboratory should have a collection of these fine quality slides for reference and teaching purposes. These PVA prepared smears are stained with Trichrome stain. Pick and choose the slides desired. **minimum of any ten (10) slides per package.**

- 4901 Entamoeba histolytica cysts and trophs
- 4907 Entamoeba hartmanni cysts and trophs
- 4909 Entamoeba coli cysts and trophs
- 4912 Endolimax nana cysts and trophs
- 4916 Dientamoeba fragilis trophs
- 4918 Giardia lamblia cysts and trophs
- 4923 Chilomastix mesnili cysts and trophs
- 4927 Trichomonas vaginalis
- 4929 Blastocystis hominis
- 4931 Iodamoeba sps. cysts and trophs
- 4945 Cryptosporidium sps. (Mod. Acid Fast Stain)
- 4953 Microsporidium sps.
- 4948 Isospora (Acid Fast)




4949 Cyclospora Stained Teaching Slides (package of 5 slides)

Helminth Ova & Larvae

These preparations of ova and larvae are excellent for teaching and for control material. Each vial contains 2ml of a formalin fixed suspension of feces containing many ova or larvae.

- 511 Trichuris sps.
- 512 Strongyloides larvae
- 515 Necator americanus
- 517 Taenia sps.
- 518 Hymenolepis nana
- 520 Schistosoma japonicum
- 521 Toxocara sps.
- 522 Ascaris lumbricoides
- 524 Schistosoma haematobium



Protozoa Suspensions & Smears

These formalin fixed suspensions and PVA unstained slides of the more commonly seen protozoa make an excellent control for every laboratory. These PVA slides have many characteristic protozoa. Each 2ml formalin vial contains many representative organisms.

PVA slides


10 SLIDE MINIMUM

- 461 Iodamoeba sps.
- 470 Entamoeba histolytica
- 471 Entamoeba hartmanni
- 472 Entamoeba coli
- 473 Endolimax nana
- 474 Dientamoeba fragilis
- 475 Giardia lamblia
- 476 Chilomastix mesnili
- 478 Blastocystis hominis



2ml Formalin Vial

- 450 Entamoeba histolytica
- 451 Entamoeba hartmanni
- 452 Entamoeba coli
- 453 Endolimax nana
- 454 Dientamoeba fragilis
- 455 Giardia lamblia
- 456 Chilomastix mesnili
- 458 Blastocystis hominis
- 460 Iodamoeba sps.
- 483 Cryptosporidium sps.
- 490 Mixed Protozoa





Acid Fast Control Slides & Stains
convenient control slides

These unique Acid Fast Control Slides simulate a clinical sputum smear with appropriate background material. Slides have two wells containing a positive and negative smear. Can be used with Kinyoun (Cat. No. 353), Ziehl-Neelsen (Kit is Cat. No. 377) or SDL's new product AcriFluor™.

- 350 Acid Fast Stain Control Slides (10)
- 351 Acid Fast Stain Control Slides (50)
- 353 AFB Stain/Kinyoun Stain Set (250ml)
- 353-A AFB Kinyoun Stain (3x250ml)
- 353-B AFB Kinyoun Decolorizer (3x250ml)
- 353-C AFB Kinyoun Counterstain (3x250ml)
- 358 Acid Fast Tissue Control Slides (10)
- 377 Ziehl-Neelson Modified AFB Kit (Stain & Brilliant Green, 250ml each)
- 384 AcriFluor™ (125ml)



Cryptosporidium or Nocardia Control Slides & Stains
many typical organisms on each slide

Cryptosporidium or Nocardia Control Slides have many typical organisms on each slide. Each package contains 10 slides, each having a positive and negative well. For the Nocardia the negative well is *Streptomyces* spp. For staining these slides SDL's product AcriFluor™ or Crypto Kinyoun Standard Stain Sets (Cat No. 373) performs well. Each Mini-Stain Set contains two 15ml dropping bottles with Stain and Counterstain and 30ml of Decolorizer. The Standard Stain Set contains three 60ml bottles.

- 370 Cryptosporidium Control Slides (10)
- 374 Nocardia Control Slides (10)
- 384 AcriFluor™ (125ml)
- 373 Cryptosporidium Kinyoun Stain Set (3x250ml)
- 373-A Cryptosporidium Kinyoun Stain (60ml)
- 373-B Cryptosporidium Kinyoun Decolorizer (60ml)
- 373-C Cryptosporidium Kinyoun Counterstain (60ml)
- 379 Cryptosporidium Kinyoun Large Stain Set (250ml)



Gram Stain Control Slides
smears on the same slide with the control organism

Our Gram Stain Control Slides are especially designed so that clinical smears can be placed on the same slide with the control organism (*Staph. aureus* and *E.coli*). Space is provided for writing data. Packaged in quantities of 10 or 50 slides.

- 360 Gram Stain Control Slides (10)
- 361 Gram Stain Control Slides (50)



Isospora Control Slides
representative of immuno-compromised patient samples

Each slide has numerous organisms for easy recognition. Staining instructions included. For staining products refer to SDL's AcriFluor or Crypto Kinyoun Standard Stain Set.

- 378 Isospora Belli Control Slides (10)



Microsporidia & Microsporidia Family Control Slides
the whole family of Microsporidia organisms

These fixed and ready to use control slides contain many representative organisms. A negative smear is also included on each slide. The Microsporidium Family Stain Control Slides bring the morphology differences of the whole family to your laboratory on one slide with a labeled sample of *Encephalitozoon hellum*, *Enterocytozoon bieneusi*, *Nosema corneum* & *Septata intestinalis*. In the package of 10, one slide is stained with Modified Trichrome Stain and the balance of unstained slides can be used as controls.

- 369 Microsporidia Control Slides (10)
- 367 Microsporidia Family Control Slides (10)



Pneumocystis Control Slides
high quality control slides

These Pneumocystis Control Slides contain many organisms in a homogenate of lung tissue similar to material received in the microbiology laboratory. Packaged in quantities of 10 slides. These slides should not be used with current monoclonal antibody fluorescent assays.

- 356 Pneumocystis Carinii Stain Control Slides (10)



Protozoa Stain Control Slides
excellent controls

These Protozoa Control Slides contain smears of a mixture of a number of the more commonly seen clinical protozoa in PVA. They are ideal for the control of Trichrome or Hematoxylin staining procedures. Now available fixed with SAF, Copper PVA, or Zinc PVA.

- 354 Protozoa Stain Control Slides (10)
- 355 Protozoa Stain Control Slides (50)
- 354-CU Protozoa Stain Control Slides fixed with Copper (10)
- 354-SAF Protozoa Stain Control Slides fixed with SAF (10)
- 354-ZC Protozoa Stain Control Slides fixed with Zinc-PVA (10)



Why invest thousands of dollars in a specialized camera microscope when you can get a ...

Microscope Video Eyepiece

affordable video eyepieces for your computer

This inexpensive color digital video camera allows the user to capture the image seen under the microscope with amazing clarity in simple .jpg format. Adapter allows camera to slide into use on most microscope eyepieces. Assembles in seconds. Line from camera fits into USB port on your computer. Easy to understand software included which lets you snap a picture and perform simple color correction. Comes with eyepiece, adapter, USB cable and software. Compatible with Windows® 95, 98, 98se, 2000 and XP. You can easily download a Macintosh® version from the software manufacturer for a small additional fee.

Captured images can be up to 640 x 480 pixels. The included software allows you to print full-sized images in various sizes and layouts and create slide shows. With an internet connection you can easily e-mail images.

Software System Requirements: Intel® Pentium® III compatible processors, Microsoft® Windows® 98SE, 2000, Me and XP, 32MB of RAM, 18MB of available hard drive space, CD-ROM drive, True Color or, HiColor display adapter and monitor; 800x600 display resolution, Windows® compatible pointing device, Microsoft® Internet Explorer 4.0 or higher



8201 Digital Video Eyepiece

ID Discs & Tablets *easy and convenient*

These Rapid ID Discs and tablets are easy to use and give easy-to-read results. Reaction times are less than 30 minutes. Weetabs are small tablets in a 10 x 75 mm tube that react with 0.5 ml of water and inoculum.

398-0485	Nitrate ID Disc	50 discs
398-0520	Oxidase ID Disc	50 discs
398-1650	Urease Weetab	28 discs
398-0915	Gram Reaction ID Disc	30 discs
398-1045	Glycosidases Weetab	28 discs
	Alkaline phosphatase	
398-1270	Glycosidases Weetab	28 discs
	Alpha-glucosidase	
398-1275	Glycosidases Weetab	28 discs
	Beta glucosidase	



398-1378	LAP ID Disc	50 discs
398-1490	Glycosidases Weetab	28 discs
	Beta galactosidase	
398-1532	PRO ID Disc	50 discs
398-1537	PYZ Weetab	28 discs
398-1538	PYR ID Disc	50 discs
398-1651	Urea/Indole Weetab	28 discs
398-9131	Butyrate ID Disc	50 discs
398-9142	Candida Screen ID Disc	25 discs
398-9145	ALN ID disc	50 discs



Humidity Chambers/Slide Moisture Chamber *convenient and affordable*

These moisture chambers are ideal for not only immunostaining but diagnostic reactions performed on slides. Each unit is divided into ten (10) compartments with a lid barrier between each slide. Chamber is fabricated out of heavy-duty polystyrene and is air tight to keep moisture in. Boxes are stackable. Comes in black or clear.



- 197BL Black Slide Moisture Chamber
- 197CR Clear Slide Moisture Chamber

McFarland Standards *priced far less than the competition*

Now you can have the same quality of McFarland Standards at a fraction of the cost. Each particle is the size of E. coli. with an extended shelf life. Each tube contains 15ml of suspension. **Expiration date of two years.**

2300	McFarland 0.5 (each)
2350	McFarland Set (one each 0.5, 1.0, 2.0, 3.0, &4.0)
2370	Colorimeter Set (one each 0.5, 10., 2.0 & 3.0)





Slide Dryer, Loop Rest, & Loop Support

helpful additions to your incinerator

These additions to your incinerator will pay for themselves many times over. The innovative Slide Dryer allows 1, 2 or 3 slides to dry quickly. The Loop Rest helps in protecting the handle of loops from overheating and melting. Loop Supports avoid burning your desk.

- 500 Slide Dryer (1 place)
- 501 Slide Dryer (2 place)
- 502 Slide Dryer (3 place)
- 700 Loop Support
- 710 Loop Rest

Slide Heater

Safely fix, dry & stain slides with a convenient heater

This low-profile slide heater is uniquely safe because it runs on a low 12 volts. Measuring 10" x 4" it is ideal to fix, dry, or stain slides. Thermostatically controlled, the temperature can be set anywhere from room temperature to 65°C. This large heater holds up to twenty (20) slides. Hot surface is made of a resin powder coated aluminum.

- 504-20 Slide Heater (20 place)
- 503 Replacement Digital Thermometer



Low Power Stain or Water Bath

This unique 12 volt high temperature bath considerably speeds up staining of larger amount of slides. It can be used for either Trichrome or Modified Trichrome procedures. It is made of easy to clean 20 gauge stainless steel. Reaches temperatures of up to 70°C with an easy-to-read and set LCD display. Comes with rack that holds up to 36 slides. Microwave resistant racks can be purchased separately. The 120V transformer that powers the unit is UL and SA approved.

- 507 Stain or Water Bath
- 9803 Microwave rack

Odor Controller *dispels odors with fresh, natural scents*

This adjustable deodorizer is a non-toxic, hypo-allergenic and flame resistant unit that smells as if the real plant was in the room. Each fragrance cartridge (made from natural plant based oils) lasts 3 months depending on room and dispenser size. The Supra unit has a light-activated auto shut-off system. Double stick foam tape is supplied for easy installation. Units do not come with odor cartridges, you must order these separately. Cartridge choices are: Apple, Crusair (orange), Lemon, Grapefruit, Hibiscus, Jasmine, Orchard Peach, Rain Forest, Sun Green, or Tropic Delite.

- 97100 Odor Controller -Tan Unit (no cartridges included)
- 97101 Odor Controller -Silver Deluxe Unit (no cartridges included)
- 97106 Odor Controller -Supra Unit (no cartridges included)
- 97110 Odor Controller Cartridges (5)

Petri Dish Desk Holders

easy to use stacking and revolving racks

Each section of the Petri Dish Desk Holder (left) and the Revolving Petri Dish Holder (right) holds 13 Petri dishes. Revolving model rotates 360°. Ideal for separation of cultures from different sources. These racks store conveniently in the corner of an incubator. Made of durable and easy to care for acrylic.

- 1814 Desk Holder (four place)
- 181 Desk Holder (seven place)
- 184 Revolving Petri Dish Holder

Teflon Bench Laminate

protects your work surfaces from stains, corrosive acids and other chemicals

This laminate protects your work surfaces from stains, corrosive acids and other chemicals. With the wipe of a cloth you have a clean bench top. Items up to 200°F can be placed on it with no surface damage. The unique protective cover is a thin film of Teflon FEP backed with peel-off sticky vinyl. Easy to apply, by cutting to the desired shape and pressing into place. Available in sizes of 15" by 25" or 9" by 29".

- 220 Teflon Bench Laminate (15" x 25")
- 221 Teflon Bench Laminate (9" x 29")

Tissue Grinders *sterile grinders great for small jobs*

These sterile disposable tissue grinders are ideal for needle biopsy specimens or small amounts of tissue meant for microbiological analysis. Grinder is made of high quality borosilicate glass with etched bottom for greater abrasion. The tube is non-cracking polypropylene and the high temperature vinyl top covers the tube and protects from breakage. Each grinder is individually packaged in a sterile, easy-to-open pouch. Tubes available in 8mm, 13mm and 17mm.

- 825-08 Tissue Grinder 8mm (12)
- 825-13 Tissue Grinder, 13mm (12)
- 825-17 Tissue Grinder, 17mm (12)
- 825-MIX Tissue Grinder (four of each size)

LAB-EL Computer Generated Labeling Program *design and print labels on your own printer*

Better identify & organize slides and test tubes when you design and print labels on your own printer. Up to 26 frequently used custom labels can be saved and LAB-EL can sequentially number. LAB-EL is an inexpensive, simple-to-use computer program that allows you to design and print on labels that perfectly fit on the end of microscope slides. You can also label aliquot-containing test tubes meant for storage. LAB-EL meets criteria of laboratory inspection by accrediting agencies. The program allows for hospital name, location and other information to be placed in the memory for easy recall. Prints a maximum of 17 letters per line, five lines per label. Consecutive numbers can be easily printed. Depending on your printer, choose from white or colored Pin-fed labels or white Laser Printer compatible label. The program is available for IBM PC or other PC compatible machines (at least 128k).



- 014 LAB-EL Program 3.5" Disk
- 005-5 Pin Feed Slide Labels-Color (5,000)
- 005-10 Pin Feed Slide Labels-Color (10,000)
- 012 White Slide Labels (5,000)
- 013 White Slide Labels (10,000)
- 018-06 Laser Slide Labels White (6,000)
- 018-12 Laser Slide Labels White (12,000)
- 018-30 Laser Slide Labels White (30,000)
- 018-50 Laser Slide Labels White (50,000)



2005 SDL Price List as of 7/1/05

SCIENTIFIC DEVICE LABORATORIES JULY 2005 PRICE LIST

Catalog #	Description	Retail
.005 -10	SLIDE LABELS-10,000 COLOR	\$46.50
.005 -5	SLIDE LABELS-5,000 COLOR	\$31.50
.012 .	WHITE SLIDE LABELS-5,000	\$28.50
.013 .	WHITE SLIDE LABELS-10,000	\$44.00
.014 .	LAB-EL PROGRAM 3.5" DISK	\$95.00
.018 -06	LASER LABELS WHITE-6,000	\$46.50
.018 -12	LASER LABELS WHITE-12,000	\$87.50
.018 -30	LASER LABELS WHITE-30,000	\$210.00
.018 -50	LASER LABELS-WHITE 50,000	\$337.50
.025 .	Atlas Parasitology CD Rom	\$120.00
.030 -1	GRAM SLIDES BX OF 72=1/2 gross *	\$26.00
.030 -2	GRAM SLIDES BX OF 72=1/2 gross *	\$26.00
.030 -3	GRAM SLIDES BX OF 72=1/2 gross *	\$26.00
.030 -999	CUSTOM GRAM SLIDE-special *	\$26.00
.041 .	CUSTOM PRINTED SLIDES /quantity 1/2 gross) *	\$26.00
.044 .	GRID SLIDES / 1/2 gross *	\$27.00
.044 -POLY	GRID SLIDES W/POLY-L-LYSINE / 1/2gross	\$32.00
.044 -SIL	GRID SLIDES W/SILANE /1/2 gross *	\$32.00
.050 -S	Sterile Slides (50 slides/case) *	\$55.00
.074 -1	LARGE AGGLUT SLIDE-SQUARE WELL (1/2 gross) *	\$40.00
.074 -2	LARGE AGGLUT SLIDES-ROUND WELL (1/2 gross) *	\$40.00
SDL100 .	QC BOARD FOR STEEL SURFACE (6)	\$42.50
SDL100 -6	QC BOARD BACKBONES ONLY-SET OF 6	\$12.50
SDL101 .	QC BOARD NON-STEEL SURFACE (6)	\$42.50
SDL121 .	ACRYLIC VIAL RACK (2 WELL/12 mm)	\$20.00
SDL122 .	ACRYLIC VIAL RACK (2 WELL/15mm)	\$20.00
SDL123 .	ACRYLIC VIAL RACK (2 WELL/17mm)	\$20.00
SDL124 .	ACRYLIC VIAL RACK (2WELL/19mm)	\$20.00
SDL125 .	ACRYLIC VIAL RACK (2 WELL/23mm)	\$20.00
SDL129 .	MULTI-TIER RACK -SET OF 2	\$23.00
SDL130 .	MYCOBACTERIA RACK-small hole	\$40.00
SDL131 .	MYCOBACTERIA RACK-large hole	\$40.00
SDL132 .	REAGENT RACK-EACH (for AP bottles)	\$40.00
SDL134 .	REAGENT RACK-each (Micros	\$40.00
SDL135 .	ACRYLIC VIAL RACK-3 WELL and hole diameter	\$26.50
SDL141 .	ACRYLIC VIAL RACK(4 WELL/12mm)	\$26.50
SDL142 .	ACRYLIC VIAL RACK(4 WELL/15mm)	\$26.50
SDL143 .	ACRYLIC VIAL RACK(4 WELL/17mm)	\$26.50
SDL144 .	ACRYLIC VIAL RACK(4 WELL/19mm)	\$26.50
SDL145 .	ACRYLIC VIAL RACK(4 WELL/23mm)	\$26.50
SDL151 .	ACRYLIC VIAL RACK (5 WELL/12mm)	\$32.50
SDL152 .	ACRYLIC VIAL RACK 5 WELL/15mm)	\$32.50
SDL153 .	ACRYLIC VIAL RACK(5 WELL/17mm)	\$32.50
SDL154 .	ACRYLIC VIAL RACK(5 WELL/19mm)	\$32.50
SDL155 .	ACRYLIC VIAL RACK(5 WELL/23mm)	\$32.50
SDL162 .	ACRYLIC VIAL RACK (6 WELL/15mm)	\$36.00
SDL165 .	ACRYLIC VIAL RACK(6 WELL/23mm)	\$36.00
SDL171 .	ACRYLIC VIAL RACK (7 WELL/12mm)	\$36.00
SDL172 .	ACRYLIC VIAL RACK(7 WELL/15mm)	\$36.00
SDL173 .	ACRYLIC VIAL RACK(7 WELL/17mm)	\$36.00
SDL174 .	ACRYLIC VIAL RACK(7 WELL/19mm)	\$36.00
SDL175 .	ACRYLIC VIAL RACK(7 WELL/23mm)	\$36.00
SDL180 .	MICROMATE-each	\$125.00
SDL181 .	PETRI DISH DESK HOLDER - 7 place	\$85.00
SDL183 .	RODAC DISH HOLDER - each	\$80.00
SDL184 .	PETRI DISH STORAGE TURNTABLE	\$77.50
SDL186 .	INCUBATOR MINI-SHELVES	\$65.00

"5" on red bckgrnd=change
 PRICES EFFECTIVE 07/01/05

7/8/2005

2005 SDL Price List as of 7/1/05

SDL188 .	Inactivating Reagent (500ml)	\$35.00
SDL189 .	Inactivating Rack -each	\$45.00
SDL190 .	TB/MYCOLOGY TUBE RACK-4 TRAY	\$135.00
SDL191 .	TB/MYCOLOGY TUBE RACK-2 TRAY	\$100.00
SDL192 .	TB Safety Rack	\$40.00
SDL193 .	SLIDE INCUBATOR	\$100.00
SDL194 .	ROOM TEMPERATURE INCUBATOR	\$140.00
SDL195 .	Desiccator/Humidifier	\$120.00
SDL196 .	FREEZER STORAGE BOX	\$45.00
SDL197 -BL	Slide Moisture Chamber (Black)	\$30.00
SDL197 -CR	Slide Moisture Chamber (clear)	\$30.00
SDL200 .	PLASTIC SLIDE INCUBATOR/ Set of 10	\$49.00
SDL200 -50	PLASTIC SLIDE INCUBATOR-SET of 50	\$176.00
SDL201 .	Rh AND TEMPERATURE DIGITAL INDICATOR-CELSIUS	\$66.00
SDL210 .	THERMOMETER-each	\$14.00
SDL211 .	THERMOMETER-SET OF 6	\$80.00
SDL217 .	HOOD VELOCITY METER	\$50.00
SDL220 .	TEFLON LAMINATE- 5 yards x25'	\$200.00
SDL221 .	TEFLON LAMINATE 9' x 29"	\$175.00
SDL250 .	TB SPLASH PROOF CONTAINER SET-250ml	\$45.00
SDL250 -999	FUNNEL & STOPPER SET REPLACEMENT (4)	\$28.50
SDL251 .	TB SPLASH PROOF CONTAINER SET-500ml	\$51.00
SDL251 -999	FUNNEL&STOPPER SET 4/REPLACEMENT	SDL28.5
SDL252 .	TB SPLASH PROOF CONTAINER SET - 1000ml	\$61.50
SDL252 -997	TB Splashproof Container No rack - 2 bottles, 2 funnels	\$40.00
SDL252 -999	FUNNEL&STOPPER SET 4/REPL	\$30.00
SDL259 .	STACKING BLOCKS-SET OF 18	SDL21.5
SDL261 .	SHRINK SEALS/400 PER CONTAINER	\$75.00
SDL262 .	SHRINK SEALS-6,000/CONTAINER	\$900.00
SDL264 .	MICROTITER SHRINK SEALS-1 gross	\$30.00
SDL265 .	SHRINK SEAL-TEST TUBE SIZE (3 gross)	\$40.00
SDL266 .	SHRINK SEAL-60MM X SDL25MM / 3 gross	\$46.50
SDL267 .	SHRINK SEALS 16 X 15 - 3 gross	\$51.75
SDL268 .	SHRINK SEALS 18 X 15 - 3 gross	\$51.75
SDL269 .	SHRINK SEALS 21 X 20 - 3 gross	\$51.75
SDL270 .	SHRINK SEALS 27 X 25 - 3 gross	\$51.75
SDL271 .	SHRINK SEALS 35 X 25 - 3 gross	\$51.75
SDL272 .	SHRINK SEAL(BDwhiskey) - 1 gross	\$51.75
SDL273 .	SHRINK SEAL (BD 8A) - 1 gross	\$51.75
SDL274 .	SHRINK SEAL(ANAEROBIC) - 1 gross	\$51.75
SDL275 .	SHRINK SEALS (DIFCO EZDRAW) - 1 gross	\$51.75
SDL276 .	SHRINK SEAL (DIFCO ESP) - 1 gross	\$51.75
SDL277 .	SHRINK SEAL(ORGANON) - 1 gross	\$51.75
SDL350 .	ACID FAST STAIN CONTROL - 10	\$19.00
SDL351 .	ACID FAST STAIN CONTROL - 50	\$83.25
SDL353 .	AFB STAIN/KINYOUN STAIN SET - 250ml	\$68.25
SDL353 -A	AFB KINYOUN STAIN -3 BOTTLES 250ml	\$94.50
SDL353 -B	AFB KINYOUN DECOLORIZER - 3 BOTTLE 250ml	\$47.25
SDL353 -C	AFB KINYOUN COUNTERSTAIN - 3 BOTTLES 250ml	\$63.00
SDL354 .	PROTOZOA STAIN CONTROL - 10	\$19.00
SDL354 -CU	PROTOZOA CONTROLS W/COPPER - 10	\$19.00
SDL354 -SAF	PROTOZOA CONTROLS W/SAF - 10	\$19.00
SDL354 -ZC	PROTOZOA CONTROL W/ZINC - 10	\$19.00
SDL355 .	PROTOZOA STAIN CONTROL - 50	\$83.25
SDL356 .	PNEUMOCYSTIS CARINII STAINED CNTRL SLIDES - 10	\$41.50
SDL357 .	YEAST CONTROL SLIDES - 50	\$75.00
SDL358 .	ACID FAST TISSUE CONTROL SLIDES - 10	\$35.00
SDL360 .	GRAM STAIN CONTROL (10)	\$19.00
SDL361 .	GRAM STAIN CONTROL (50)	\$83.25
SDL362 .	SPORE STAIN CONTROL (10)	\$19.00

2005 SDL Price List as of 7/1/05

SDL363 .	SPORE STAIN CONTROL (50)	\$83.25
SDL367 .	MICROSPORIDIUM FAMILY CONTROL SLIDES - 10	\$52.00
SDL369 .	MICROSPORIDIA CONTROL SLIDE - 10	\$33.75
SDL370 .	CRYPTOSPORIDIUM CONTROL SLIDE -10	\$33.75
SDL373 .	CRYPTOSPORIDIUM KINYOUN STAIN SET - 60ml	\$52.00
SDL373 -A	CRYPTOSPORIDIUM KINYOUN STAIN - 60ml	\$19.50
SDL373 -B	CRYPTOSPORIDIUM KINYOUN DECOLORIZER- 60ml	\$19.50
SDL373 -C	CRYPTOSPORIDIUM KINYOUN COUNTERSTAIN - 60ml	\$19.50
SDL374 .	NOCARDIA/STREPTOMYCES SLIDE -10	\$36.50
SDL377 .	ZIEHL-NEELSON MODIFIED AFB KIT	\$68.25
SDL378 .	ISOSPORA BELLI-10 SLIDES	\$40.00
SDL379 .	CRYPTO KINYOUN LARGE STAIN SET - 250ml	\$68.25
SDL379 -3	STAIN RACK FOR 250ML BOTTLE-each	\$25.00
SDL379 -B	Decolorizer - 6 pack 250ml each	\$55.00
SDL381 .	ZIEHL NEELSON BRILLIANT GREEN -250ml	\$30.00
SDL384 .	ACRI-FLUOR STAIN KIT-125ml	\$48.00
SDL390 .	HEMO P Stain Set (40ml)	\$38.50
SDL390 -1	HEMO P Reagent A Replacment (40ml)	\$22.00
SDL390 -2	HEMO P Reagent B Replacment (40ml)	\$22.00
SDL390 -3	HEMO P Fixative Replacment (40ml)	\$12.00
SDL391 .	HEMO P Stain Set (250ml)	\$93.00
SDL391 -1	HEMO P Reagent A Replacment (250ml)	\$55.00
SDL391 -2	HEMO P Reagent B Replacment (250ml)	\$55.00
SDL391 -3	HEMO P Fixative Replacment (250ml)	\$27.50
SDL391 -F	HEMO P Fecal Stain Set	\$93.50
SDL392 -250	Gram Stain Set (250ml)	\$45.00
SDL392 -60	Gentian Violet per 6 btls of 250ml each	\$72.00
SDL392 -70	Gram's Iodine per 6 btls of 250ml each	\$45.00
SDL392 -80	Gram's Decolorizer 6 btls 250ml each	\$40.00
SDL392 -90	Gram's Safranin 6 btls of 250ml each	\$72.00
SDL398 -0485	Nitrate ID Disc (50 discs)	\$14.00
SDL398 -0520	Oxidase ID Disc (50 discs)	\$11.00
SDL398 -0915	Gram Reaction ID Disc (30 discs)	\$23.00
SDL398 -1045	Glycosidases Weetab Alkaline phosphatase(28 discs)	\$14.00
SDL398 -1270	Glycosidases Weetab Alpha-glucosidase(28 discs)	\$14.00
SDL398 -1275	Glycosidases Weetab Beta glucosidase(28 discs)	\$14.00
SDL398 -1378	LAP ID Disc (50 discs)	\$23.00
SDL398 -1490	Glycosidases Weetab Beta galactosidase (28 discs)	\$14.00
SDL398 -1532	PRO ID Disc (50 discs)	\$23.00
SDL398 -1537	PYZ Weetab (28 discs)	\$14.00
SDL398 -1538	PYR ID Disc (50 discs)	\$28.00
SDL398 -1650	Urase (28 discs)	\$11.00
SDL398 -1651	Urea/Indole Weetab (28 discs)	\$11.00
SDL398 -9131	Butyrate ID Disc (50 discs)	\$39.00
SDL398 -9142	Candida Screen ID Disc (25 discs)	\$21.50
SDL398 -9145	ALN ID disc (50 discs)	\$23.00
SDL430 .	GIEMSA/WRIGHT CONTROL -10	\$41.50
SDL431 .	IRON HEMATOLOGY STAIN CONTROL -10	\$41.50
SDL432 .	LEUCOCYTE ALKALINE PHOSPH	\$37.50
SDL433 .	IRON CONTROL SLIDES -50	\$140.00
SDL448 .	POT LUCK PARASITLOGY SPECIMENS	\$25.00
SDL450 .	E. HISTOLYTICA-2CC VIAL	\$50.00
SDL451 .	E. HARTMANNI-2CC VIAL	\$50.00
SDL452 .	E. COLI-2CC VIAL	\$50.00
SDL453 .	E. NANA-2CC VIAL	\$50.00
SDL454 .	D. FRAGILIS-2CC VIAL	\$50.00
SDL455 .	GIARDIA LAMBLIA-2CC VIAL	\$50.00
SDL456 .	C. MESNILI-2CC VIAL	\$50.00
SDL458 .	BLASTOCYSTIS HOMINIS-2CC VIAL	\$50.00
SDL460 .	IODAMOEBIA SPS-2CC VIAL	\$50.00

2005 SDL Price List as of 7/1/05

SDL461 .	IODAMOEBIA SPS-PVA UNSTAINED SLIDE	\$50.00
SDL469 .	DIROFILARIA IMMITIS (5)	\$45.00
SDL470 .	E. HISTOLYTICA-PVA UNSTND - 10	\$50.00
SDL471 .	E. HARTMANNI-PVA UNSTND - 10	\$50.00
SDL472 .	E. COLI-PVA UNSTND (10)	\$50.00
SDL473 .	E. NANA-PVA UNSTND (10)	\$50.00
SDL474 .	D. FRAGILIS-PVA UNSTND - 10	\$50.00
SDL475 .	GIARDIA LEMBLIA-PVA UNSTAINED SLIDE - 10	\$50.00
SDL476 .	CHILOMASTIX MESNILI-PVA UNSTAINED - 10	\$50.00
SDL478 .	BLASTOCYSTIS HOMIS-PVA UNSTAINED - 10	\$50.00
SDL483 .	CRYPTOSPORIDIUM SPS-2CC VIAL	\$50.00
SDL490 .	MIXED PROTOZOA-2CC VIAL	\$50.00
SDL491 .	PLASMODIUM VIVAX-BLOOD SMEAR -5 SLIDES	\$45.00
SDL492 .	PLASMODIUM MALARIAE-BLOOD - 5 SLIDES	\$45.00
SDL493 .	PLASMODIUM FALCIPARUM-BLOOD - 5 SLIDES	\$45.00
SDL494 .	TRYPANOSOMA RHODESIENSE - 5 SLIDES	\$45.00
SDL496 .	LOA-LOA 5 SLIDES	\$45.00
SDL497 .	BRUGIA MALAYI-5 SLIDES	\$45.00
SDL500 .	SLIDE DRYER-1 PLACE	\$35.00
SDL501 .	SLIDE DRYER-2 PLACE	\$40.00
SDL502 .	SLIDE DRYER-3 PLACE	\$55.00
SDL504 -16	SLIDE HEATER W/ TRANSFORMER-high heat	\$300.00
SDL504 -20	SLIDE HEATER W/ TRANSFORMER	\$265.00
SDL510 .	E. VERMICULARIS-2CC VIAL	\$60.00
SDL511 .	TRICHURA TRICHURIS-2CC VIAL	\$60.00
SDL512 .	STRONGYLOIDES LARVAE-2CC VIAL	\$60.00
SDL513 .	SHISTOSOMA MANSONI-2CC VIAL	\$60.00
SDL515 .	NECATOR AMERICANUS-2CC VIAL	\$60.00
SDL517 .	TAENIA SPS-2CC VIAL	\$60.00
SDL518 .	HYMENOLOPIS NANA -2CC VIAL	\$60.00
SDL519 .	D. LATUM- 2CC VIAL	\$60.00
SDL520 .	SHISTOSOMA JAPONICUM-2CC VIAL	\$60.00
SDL521 .	TOXOCARA SPS-2CC VIAL	\$60.00
SDL522 .	ASCARIS LUMBRICOIDES-2CC VIAL	\$60.00
SDL524 .	S. HAEMATOBIIUM-2CC VIAL	\$60.00
SDL526 .	FASCIOLA HEPATICA-2CC VIAL	\$60.00
SDL550 .	TRICHINELLA SPIRALIS-MUSCLE - 5 SLIDES	\$45.00
SDL552 .	LEISHMANIA DONOVANI - 5 slides	\$45.00
SDL553 .	SHISTOSOMA MANSONI-5 SLIDES	\$45.00
SDL560 .	ENTAMOEBIA HISTOLOYTICA-CECUM - 5 SLIDES	\$45.00
SDL581 .	EXPLOSTION-PROOF STOOL CONTAINERS SET OF 100	\$72.00
SDL582 .	PINWORM COLLECTOR 25 vial	\$40.00
SDL614 -01	NOTOX- 1 GALLON	\$24.50
SDL614 -05	NOTOX 5 GALLON CUBE	\$105.00
SDL614 -120	NOTOX 120ML CONTAINER 50/case	\$34.00
SDL614 -15	NOTOX 15 ML CONTAINERS 50/case	\$29.00
SDL614 -180	NoTox 180ml 24/box	\$20.00
SDL614 -240	NoTox 240ml box of 24	\$26.00
SDL614 -250	NoTox 250ml - box of 24	\$26.00
SDL614 -30	NOTOX 30 ML CONTAINERS 50	\$24.50
SDL614 -60	NOTOX 60 ML CONTAINERS 50	\$31.00
SDL614 -90	NOTOX 90 ML CONTAINERS 5	\$32.50
SDL615 -1	PARACLEAR - 1 GALLON UNIT	\$35.00
SDL615 -4	PARACLEAR 4 X 1 GALLON UNIT	\$125.00
SDL615 -5	PARACLEAR 5 GALLON CUBE	\$130.00
SDL616 -100	SUREMOUNT - 100 ML BOTTLE	\$12.50
SDL616 -400	SUREMOUNT 4 X 100 ML BOTTLE	\$45.00
SDL616 -500	SUREMOUNT 500 ML BOTTLE	\$45.00
SDL617 .	DETOX NEUTRALIZER 10 LB JAR	\$50.00
SDL617 -01	DETOX NEUTRALIZER/BOX OF 8 CONTAINERS	\$32.50

"5" on red bckgrnd=change
PRICES EFFECTIVE 07/01/05

7/8/2005

2005 SDL Price List as of 7/1/05

SDL617 -011	DETOX NEUTRALIZER/BX OF 4 CONTAINERS	\$37.50
SDL617 -04	DETOX NEUTRALIZER 4 x 10LB JARS	\$180.00
SDL617 -100	DeToX NEUTRALIZER (100 VIALS)	\$45.00
SDL617 -99	DETOX CONFIRMATORY TEST - 8 tests	\$30.00
SDL618 -100	Floromount - 100ml bottle	\$15.00
SDL618 -400	Floromount - 4 x 100ml bottle	\$45.00
SDL618 -500	Floromount - 500ml bottle	\$50.00
SDL619 .	ENHANCING FLUID (250 ML)	\$10.00
SDL619 -6	ENHANCING FLUID - 6 BOTTLES - 250ml	\$50.00
SDL641 .	FECAL CONCENTRATOR KIT-50	\$62.50
SDL641 -P	FECAL CONCENTRATOR KIT FOR PARAWASH	\$49.00
SDL642 .	FINE MESH FILTER FUNNELS - 250	\$170.00
SDL643 .	CENTRIFUGE TUBES WITH CAP - 250	\$36.00
SDL644 .	MUCOLYTIC DIGESTING AGENT - 15ml	\$20.50
SDL645 .	FECAL CONCENTRATOR RACK - 6 place rack	\$35.00
SDL660 .	SNAP N' STAIN IRON STAIN SET (12)	\$35.00
SDL660 -9	COUNTERSTAIN ONLY	\$15.00
SDL661 .	SNAP N' STAIN IRON - 50 UNITS	\$67.50
SDL662 .	NUCLEAR FAST RED COUNTERSTAIN - 250ML	\$25.00
SDL667 .	SNAP N' DIGEST 75ML - 10 BOTTLES	\$48.00
SDL668 .	SNAP N'DIGEST 150ML - 10 BOTTLES	\$80.00
SDL669 -100	TB BUFFER 100ML /10 BOTTLES	\$165.00
SDL669 -500	TB BUFFER 500ML /1 BOTTLE	\$65.00
SDL670 -02	STREPTOMYCES GRISEUS-PER SLIDE	\$10.00
SDL670 -03	NOCARDIA- PER SLIDE	\$10.00
SDL671 -01	BLASTOSCHIZOMYCES CAPITATUS-PER SLIDE*	\$10.00
SDL671 -02	CANDIDA ALBICANS-PER SLIDE	\$10.00
SDL671 -03	CRYPTOCOCCUS NEOFORMANS-PER SLIDE	\$10.00
SDL671 -04	GEOTRICHUM CANDIDUM-PER SLIDE	\$10.00
SDL671 -05	HAMSENULA-PER SLIDE	\$10.00
SDL671 -07	SPOROTHRIX SCHENCKII-PER SLIDE	\$10.00
SDL672 -04	MUCOR SPS-PER SLIDE	\$10.00
SDL672 -06	RHIZOPUS-PER SLIDE	\$10.00
SDL672 -07	RHIZOMUCOR-PER SLIDE	\$10.00
SDL672 -09	SYNCEPHALASTRUM-PER SLIDE	\$10.00
SDL673 -01	ALTERNARIA-PER SLIDE	\$10.00
SDL673 -02	BIPOLARIS-PER SLIDE	\$10.00
SDL673 -03	CHAETOMIUM-PER SLIDE	\$10.00
SDL673 -04	CLADOSPORIUM-PER SLIDE	\$10.00
SDL673 -05	CURVULARIA-PER SLIDE	\$10.00
SDL673 -06	DRECHSLERA-PER SLIDE	\$10.00
SDL673 -10	NIGROSPORA -PER SLIDE	\$10.00
SDL673 -12	PHILALOPHORA VERRUCOSA-PER SLIDE	\$10.00
SDL673 -15	SCEDOSPORIUM APOSPERMUM-PER SLIDE	\$10.00
SDL673 -16	SCEDOSPORIUM PROLIFICANS-PER SLIDE	\$10.00
SDL673 -18	ULOCLADIUM-PER SLIDE	\$10.00
SDL673 -19	WANGIELLA DERMITITIDIS-PER SLIDE	\$10.00
SDL674 -01	EPIDERMOPHYTON FLOCCOSUM-PER SLIDE	\$10.00
SDL674 -02	MICROSPORIUM AUDOUINII--PER SLIDE	\$10.00
SDL674 -03	M. CANIS -PER SLIDE	\$10.00
SDL674 -07	M. GYPSEUM-PER SLIDE	\$10.00
SDL674 -08	M. VANBREUSEGHEMII-PER SLIDE	\$10.00
SDL674 -09	T. AJELLOI-PER SLIDE	\$10.00
SDL674 -11	T. EQUINUM-PER SLIDE	\$10.00
SDL674 -12	T. MEGNINII -PER SLIDE	\$10.00
SDL674 -13	T. MENTAGROPHYTES-PER SLIDE	\$10.00
SDL674 -14	T. RUBRUM-PER SLIDE	\$10.00
SDL674 -15	T. SCHOENLEINI-PER SLIDE	\$10.00
SDL674 -16	T. TERRESTRE-PER SLIDE	\$10.00
SDL674 -17	T. TONSURANS-PER SLIDE	\$10.00

"5" on red bckgrnd=change
PRICES EFFECTIVE 07/01/05

7/8/2005

2005 SDL Price List as of 7/1/05

SDL675 -01	ACREMONIUM SP.-PER SLIDE	\$10.00
SDL675 -06	ASPERGILLUS FLAVUS-PER SLIDE	\$10.00
SDL675 -07	ASPERGILLUS FURNIGATIS-PER SLIDE	\$10.00
SDL675 -09	ASPERGILLUS NIGER-PER SLIDE	\$10.00
SDL675 -16	CHYRSOSPORIUM-PER SLIDE	\$10.00
SDL675 -18	FUSARIUM-PER SLIDE	\$10.00
SDL675 -19	PAECILOMYCES-PER SLIDE	\$10.00
SDL675 -20	PENICILLIUM-PER SLIDE	\$10.00
SDL675 -32	SCOPULARIOPSIS-PER SLIDE	\$10.00
SDL675 -33	SEPEDONIUM SP-PER SLIDE	\$10.00
SDL685 .	CURVED NEEDLE WOODEN HANDLE-12/PKG	\$20.00
SDL686 .	STRAIGHT NEEDLE WOODEN HANDLE- 12/PKG	\$20.00
SDL687 .	ALUMINIUM MYCOLOGY NEEDLES	\$30.00
SDL700 .	LOOP SUPPORT - each	\$25.00
SDL710 .	LOOP REST - each	\$28.00
SDL714 .	ANAEROBIC INCUBATOR SMALL STAINLESS STEEL JAR	\$60.00
SDL715 .	ANAEROBIC INCUBATOR LARGE STAINLESS STEEL JAR	\$70.00
SDL716 .	ANAEROBIC INCUBATOR SMALL STAINLESS STEEL JAR	\$70.00
SDL717 .	ANAEROBIC INCUBATOR LARGE STAINLESS STEEL JAR	\$80.00
SDL720 -10	RACK FOR ANAEROBIC JAR - 9 PLATE - each	\$20.00
SDL722 -S	REPLACEMENT SEALS FOR STAINLESS(SET OF 5)	\$18.50
SDL724 -14	RACK ANAEROBIC JAR - 12 PLATE each	\$22.00
SDL725 .	ANAEROGEN - BOX 10	\$34.50
SDL726 .	ANAEROGEN INDICATOR STRIP - 100 STRIPS	\$46.00
SDL731 .	ZINC SULFATE SOLUTION-each Specific Gravity 1.18	\$20.00
SDL732 .	ZINC SULFATE SOLUTION(6) Specific Gravity 1.18	\$100.00
SDL733 .	ZINC SULFATE SOLUTION-each Specific Gravity 1.20	\$20.00
SDL734 .	ZINC SULFATE SOLUT(6) Specific Gravity 1.20	\$100.00
SDL735 .	PVA SOLUTION - 500 ML	\$33.00
SDL736 .	PVA SOLUTION-6 BOTTLES 500ml	\$175.00
SDL740 .	MICROSCOPE CARE KIT	\$75.00
SDL741 .	LENS CLEANER (PKG 6 x 6.5ml)	\$24.50
SDL742 .	IMMERSION OIL (PKG 6 x 6.5ml)	\$42.50
SDL743 .	IMMERSION OIL SOLVENT(6 X 6.5ml)	\$24.50
SDL744 .	LENS PAPER (PKG 6)	\$19.50
SDL745 .	FUNGI-TAPE (ONE ROLL)	\$27.50
SDL745 -1	MYCOLOGY MOUNTING KIT-BLUE	\$52.50
SDL745 -9	MYCOLOGY MOUNTING KIT-RED	\$52.50
SDL746 .	MINI GRIP-IT (EACH)	\$25.00
SDL746 -03	MINI GRIP-IT RACK (HOLDS 3)	\$15.00
SDL746 -73	MINI GRIP-IT + RACK /ONE OF EACH	\$35.00
SDL748 .	GRIP-IT (SET OF TWO)	\$30.00
SDL750 .	MYCOLOGY STARTER KIT-BLUE	\$78.00
SDL750 -9	MYCOLOGY STARTER KIT-RED	\$78.00
SDL751 .	10% KOH (PKG 6x3.5 dropping bottles)	\$26.50
SDL752 .	MYCOPERM RED - 12 ML	\$30.00
SDL753 .	MYCOPERM-BLUE 12 ML/BOTTLE	\$30.00
SDL754 .	INDIA INK- 6 x 3.5ml bottles	\$45.00
SDL755 .	India Ink Capsule Postitive Control with 6x 3.5ml india ink	\$70.00
SDL756 -1	India Ink Negative Control	\$15.00
SDL756 .	India Ink Capsule Positive Control	\$30.00
SDL759 .	Bright-Fluor (6 x 3.5ml bottle)	\$35.00
SDL761	Trehalose Slides (10)	\$60.00
SDL772 .	METHYLENE BLUE - 6 x 3.5ml BOTTLES	\$25.00
SDL782 .	TRICHROME STAIN (250 ML)	\$35.00
SDL783 .	TRICHROME STAIN (6 x 250ml)	\$200.00
SDL785 .	SNAP N' STAIN FOR PARASITE WET MOUNTS (6)	\$32.50
SDL790 .	TRICHROME STAINING STATION	\$127.50
SDL792 .	IODINE REPLACEMENTS-SET OF 20	\$25.00
SDL795 .	TRICHROME STAINQUICK SYSTEM	\$170.00

"5" on red bckgrnd=change
PRICES EFFECTIVE 07/01/05

7/8/2005

2005 SDL Price List as of 7/1/05

SDL796 .	STAINQUICK TRICHROME STAIN &DECOLOR -250ml ea	\$56.00
SDL820	Sampling Grids Non -Sterile (50)	\$85.00
SDL821	Sampling Grids Sterile individually wrapped (50)	\$165.00
SDL825 -08	Tissue Grinder 8mm sterile	\$60.00
SDL825 -13	Tissue Grinder 13mm sterile	\$60.00
SDL825 -17	Tissue Grinder 17mm sterile	\$60.00
SDL825 -MIX	Tissue Grinder 4 each of 3 sizes sterile	\$60.00
SDL902 .	SLIDE FORCEPS (PKG 10)	\$16.00
SDL902 -100	SLIDE FORCEPS- PKG OF 100	\$145.00
SDL902 -50	SLIDE FORCEPS-PKG OF 50	\$75.00
SDL903 .	STERILE INSTRUMENT PACKS (12)	\$59.75
SDL903 -999	STERILE SCISSORS (each)	\$2.75
SDL903 R12	Sterile Instrument Pack with Retractable Scalpel (12)	\$85.00
SDL904 .	STERILE FORCEPS (25/LOT)	\$45.00
SDL904 NS	NON STERILE FORCEPS (100)	\$100.00
SDL905 .	CEL-TEK ACID FAST CONTRL (10)	\$31.00
SDL906 .	STERILE INSTRUMENTS W/BOARD (10)	\$74.50
SDL915 .	CEL TEK AMYLOID SLIDES-10	\$31.00
SDL925 .	CEL-TEK GRAM STAIN CONTRO	\$31.00
SDL930 .	CEL-TEK IRON STAIN CONTRO	\$31.00
SDL931 .	IRON NEGATIVE CONTROLS-TI	\$31.00
SDL940 .	CEL-TEK MUCIN STAIN CONTR	\$31.00
SDL941 .	CEL-TEK MELANIN STAIN CON	\$31.00
SDL950 .	CEL-TEK PAS STAIN CONTROL	\$31.00
SDL951 .	CEL-TEK DIASTASE DIGESTION CONTRO	\$31.00
SDL960 .	CEL-TEK METHENAMINE SILVER-FUNGUS	\$31.00
SDL967 .	COMBO REAGENT RACK 7.5 &	\$60.00
SDL968 .	ADD ON RACK-six 3.5 ml holes	\$15.00
SDL969 .	REAGENT RACK	\$65.00
SDL970 .	CEL-TEK RETICULUM CONTROL	\$31.00
SDL971 .	MINI IMMUNO RACK	\$40.00
SDL972 .	IMMUNO RACK	\$65.00
SDL973 .	PARAFFIN TRAY FOR MODEL 820 - each	\$45.00
SDL974 .	PARAFFIN TRAY 2030/2040 - each	\$45.00
SDL975 .	PARAFFIN TRAY FOR NEW A/O each	\$45.00
SDL976 .	PARAFFIN TRAY FOR OLD A/O	\$45.00
SDL979 .	REAGENT RACK -40 17ML HOL	\$65.00
SDL980 .	CEL-TEK TRICHROME STAIN SET	\$31.00
SDL983 .	8oz. Specimen Container - 100	\$65.00
SDL985 .	16 oz. Specimen Container - 100	\$65.00
SDL987 .	32 oz. Specimen Container - 100	\$73.50
SDL989 .	64 oz. Specimen Container - 50	\$77.00
SDL993 .	WATER DISTILLER (1 GALLON)	\$275.00
SDL993 -9	CARBON REPLACEMENTS FOR WATER DISTILLER (6)	\$72.50
SDL1814 .	PETRI DISH DESK HOLDER - 4 place	\$65.00
SDL1917 .	STAINING RACK & CONTAINER	\$100.00
SDL1918 .	CONTAINER REPLACEMENTS-SET OF 6	\$45.00
SDL1951 .	STAINQUICK STAIN SYSTEM - MICROSPORIDIA	\$200.00
SDL1955 .	STAINQUICK MYCOBACTERIA FOR TISSUE	\$125.00
SDL2300 .	McFarland 0.5	\$12.50
SDL2350 .	McFarland Set	\$60.00
SDL2370 .	Colorimeter Set	\$48.00
SDL2584 .	LARGE STAINING JAR -EACH	\$27.00
SDL2585 .	SMALL STAINING JAR - EACH	\$16.50
SDL2586 .	5 ML VIALS W/CAPS EMPTY - 1000	\$51.00
SDL2587 .	20ML VIALS W/CAPS EMPTY-1000	\$98.00
SDL2588 .	30ML VIALS W/CAPS EMPTY-1000	\$112.00
SDL2589 .	CASSETTE FILE DRAWER SET	\$30.00
SDL3025 .	AUTOCLAVE POUCH 3 1/2" x 5 1/4" - 100/box	\$18.50
SDL3026 .	AUTOCLAVE POUCH(100) 3 1/2" x 8 1/2" - 100/box	\$20.00

"5" on red bckgrnd=change
PRICES EFFECTIVE 07/01/05

7/8/2005

2005 SDL Price List as of 7/1/05

SDL3028 .	AUTOCLAVE POUCH 5 1/4" x 10 1/4" - 100/box	\$26.50
SDL3029 .	AUTOCLAVE POUCH 7 1/2" x 13 1/2" - 100/box	\$33.00
SDL3030 .	AUTOCALVE POUCH 12" x 15 1/2" - 100/box	\$51.00
S3800 -010	Kwik-Pour Agar, Bacteriological (12 tubes)	\$24.50
S3800 -020	Kwik-Pour Brilliant Green Agar (12 tubes)	\$24.50
S3800 -030	Kwik-Pour Brain Heart Infusion Agar (12 tubes)	\$24.50
S3800 -050	Kwik-Pour Desoxycholate Lactose Sucrose Agar (12)	\$24.50
S3800 -060	Kwik-Pour ENDO Agar (12 tubes)	\$24.50
S3800 -070	Kwik-Pour Hektoen Enteric Agar (12 tubes)	\$24.50
S3800 -080	Kwik-Pour Levine EMB Agar (12 tubes)	\$24.50
S3800 -090	Kwik-Pour Mannitol Salt Agar (12 tubes)	\$24.50
S3800 -100	Kwik-Pour Nutrient Agar (12 tubes)	\$24.50
S3800 -110	Kwik-Pour Plate Count Agar (12 tubes)	\$24.50
S3800 -120	Kwik-Pour Potato Dextrose Agar (12 tubes)	\$24.50
S3800 -130	Kwik-Pour Potato Flake Agar (12 tubes)	\$24.50
S3800 -140	Kwik-Pour Sabourauds Dextrose Agar (12 tubes)	\$24.50
S3800 -150	Kwik-Pour Sorbitol MacConkeys (12 tubes)	\$24.50
S3800 -160	Kwik-Pour SS Agar (12 tubes)	\$24.50
S3800 -170	Kwik-Pour Standard Plate Count Agar (12 tubes)	\$24.50
S3800 -180	Kwik-Pour Tryptone Soy Agar (12 tubes)	\$24.50
SDL4120 .	PUS CELLS HEMATOXYLIN - each	\$10.00
SDL4449 .	STAINED SLIDES - each	\$10.00
SDL4901 .	ENTAMOEBA HISTOLYTICA SLIDE cysts & trophs	\$6.00
SDL4907 .	ENTAMOEBA HARTMANII SLIDE	\$6.00
SDL4909 .	ENTAMOEBA COLI SLIDE	\$6.00
SDL4912 .	ENDOLIMAX NANA SLIDE	\$6.00
SDL4916 .	DIENTAMOEBA FRAGILIS SLIDE	\$6.00
SDL4918 .	GIARDIA LAMBLIA SLIDE	\$6.00
SDL4923 .	CHILOMASTIX MESNILI SLIDE	\$6.00
SDL4925 .	P. MALARIAE-GIEMSA THICK - 5 SLIDES	\$45.00
SDL4927 .	TRICHOMONAS VAGINALIS SLI	\$6.00
SDL4929 .	BLASTOCYSTIS HOMINIS SLID	\$6.00
SDL4930 .	P. FALCIPARUM (THICK DROP) - 5 SLIDES	\$45.00
SDL4931 .	IODAMOEBA SPS SLIDE	\$6.00
SDL4936 .	BABESIA SPS.- 5 SLIDES	\$45.00
SDL4939 .	ACAPLASMA SPS SLIDES (5)	\$45.00
SDL4945 .	CRYPTOSPORIDIUM SPS SLIDE	\$6.00
SDL4948 .	ISOSPORA BELLI SLIDE	\$6.00
SDL4949 .	CYCLOSPORA SLIDES-STAINED - 5 SLIDES	\$47.00
SDL4950 .	CYCLOSPORA STAIN QUICK KIT - 250ml	\$150.00
SDL4953 .	MICROSPORIDIUM STAINED SLIDE	\$6.00
SDL6000 .	MYCOEASE	\$90.00
SDL6014 .	BIOHAZARD WORK KIT /SET of 12	\$160.00
SDL6016 .	DISPOSABLE CUTTING BOARDS (24)	\$50.00
SDL6020 .	Kwik-Smear Slide for KOH prep (15)	\$25.00
SDL6025 .	KOH Fungal Control, 1cc vial	\$25.00
SDL6030 .	PARASAFE -24 VIALS	\$35.00
SDL6031 .	PARASAFE & CLEAN VIAL SET - 12 sets	\$28.00
SDL6032 .	SLIDE FIXATIVE - 6.5 ML	\$30.00
SDL6036 .	CAREY BLAIR ONLY/24 VIAL	\$22.50
SDL6037 .	PARASAFE & CAREY BLAIR SET - 12 sets	\$29.00
SDL6038 .	CLEAN VIALS - SET OF 24 vials	\$9.00
SDL6040 -100	MINI PARASAFE FIXATIVE-100 vials	\$135.00
SDL6046 -100	MINI PARASAFE W/CONCENTRATORS - 100	\$170.00
SDL6046 -25	MINI PARASAFE W/CONCENTRATORS - 25	\$48.00
SDL6051 .	TUBES AND CAPS ONLY 25/BOX	\$10.00
SDL6052 .	FECAL PARASITE CONCENTRATORS 15ml - 100	\$98.00
SDL6056 .	FECAL PARSITE CONCENTRATOR 50ml - 100	\$162.50
SDL6060 .	SAF Vials (24)	\$33.50
SDL6065 -100	SAF Mini Vials (100)	\$110.00

"5" on red bckgrnd=change
PRICES EFFECTIVE 07/01/05

7/8/2005

2005 SDL Price List as of 7/1/05

SDL6121 .	PARAWASH - 6 BOTTLES 250ml	\$37.50
SDL6671 .	BUFFER PACKETS-5 PACKETS	\$20.00
SDL6790 .	Mycology CD-Rom	\$200.00
SDL7134 .	PARASITOLOGY 2" X 2" SET A CD-ROM	\$110.00
SDL7135 .	PARASITOLOGY 2" X 2" SET B CD-ROM	\$95.00
SDL7136 .	PARASITOLOGY 2" X 2" SET C CD-ROM	\$85.00
SDL7137 .	PARASITOLOGY 2" X 2" SET A+B+C CD-ROM	\$250.00
SDL7140 .	THERMOTOTE12 ADJUSTABLE WITH TEMP PORT	\$220.00
SDL7141 .	THERMOTOTE24 ADJUSTABLE WITH TEMP PORT	\$240.00
SDL7145 .	THERMOTOTE12 ADJUSTABLE WITH DISPLAY	\$265.00
SDL7146 .	THERMOTOTE24 ADJUSTABLE WITH DISPLAY	\$285.00
SDL7148 .	ThermoTote12 adjustable w/display and BATTERY	\$620.00
SDL7149 .	ThermoTote12 w/display and BATTERY	\$600.00
SDL7151 .	DELUXE THERMOTOTE INCUBATOR	\$440.00
SDL7152 -01	multiple female plug	\$55.00
SDL7153 .	Calibration Meter	\$135.00
SDL7154 .	BATTERY PACK 12 V WITH CHARGER	\$400.00
SDL7155 INTL	220v/240v converter international	\$90.00
SDL7155 USA	AC/DC CONVERTER 120V - EACH for stdn U.S. plugs	\$65.00
SDL7156 .	BATTERY GUARD FOR THERMOTOTE	\$55.00
SDL7157 -12	ThermoTote12 Petri Dish Rack (4 RACKS)	\$75.00
SDL7157 -24	ThermoTote24 Petri Dish Rack (4 RACKS)	\$75.00
SDL7158 .	SHELVES FOR THERMOTOTE (2)	\$25.00
SDL7159 .	PLASTIC TEST TUBE RACKS (4)	\$100.00
SDL7160 .	CIGARETTE LIGHTER ADAPTER	\$30.00
SDL7162 .	CIGARETTE LIGHTER ADAPTER for small & med Thermototes	\$30.00
SDL7840 .	MICROSPORIDIA KIT-MOD TRICHROME - 250ml	\$100.00
SDL7841 .	MICROSPORIDIA STAIN-250ML	\$55.00
SDL7842 C	CYCLOSPORA DECOLORIZER 250ml	\$17.50
SDL7842 M	MICROSPORIDIA DECOLORIZER - 250ml	\$20.00
SDL7843 .	MODIFIED TRICHROME STAIN KIT-250ml	\$70.00
SDL8200	Video Eyepeice and cable	\$220.00
SDL8201 .	Digital Video EyePiece, Cable & Program	\$220.00
SDL8202	Video Eyepeice, cable & Monitor	\$375.00
SDL8905 .	MYCOBACTERIUM TUBERCULOSIS - 5 slides	\$45.00
SDL8925 .	GRAM POSITIVE & NEG BACTERIA - 5 slides	\$45.00
SDL8960 .	FUNGUS (CANDIDA ALBICANS) - 5 slides	\$45.00
SDL8967 .	TOXOPLASMA GONDI-5 SLIDES	\$45.00
SDL8969 .	ACATHAMOEBIA SPS. - 5 slides	\$45.00
SDL9001 .	DISSECTING SCISSOR 5 1/2" (6)	\$27.50
SDL9002 .	POINTED SCISSORS 4 1/2" (6)	\$19.00
SDL9003 .	STERILE SCAPEL #11 BLADE- 100.pack	\$147.00
SDL9003 R100	Sterile Safety Scalpel with Retractable Blade (box of 100)	\$300.00
SDL9004 .	FINE TIPPED FORCEPS-PKG of 12	\$20.00
SDL9004 -100	FINE TIPPED FORCEPS-PKG of 100	\$125.00
SDL9004 -50	FINE TIPPED FORCEPS-PKG of 50	\$70.00
SDL9004 -S	STERILE FINE TIPPED FORCEPS - 100	\$260.00
SDL9007 .	STRAIGHT FORCEPS 4 1/2" (6)	\$19.00
SDL9008 .	STRAIGHT FORCEPS 5 1/2" (6)	\$19.00
SDL9010 .	BLUNT FORCEPS 5 1/2" - each	\$20.00
SDL9012 .	POINTED FORCEPS 6 1/2" - each	\$20.00
SDL9050 .	DISSECTION BOARD - EACH	\$30.00
SDL9050 -999	RULER AND CLIP REPLACEMENT	\$8.50
SDL9711 .	MINI STAINING STATION-2 PLACE	\$30.00
SDL9712 .	MINI STAINING STATION-6 PLACE	\$55.00
SDL9713 .	STAINING RACK W/10 CONTAINERS	\$93.00
SDL9713 -999	MINI-STAINING STATION-12 PLACE	\$120.00
SDL9714 .	STAINING CONTAINERS SET of 6	\$49.50
SDL9750 -20	Fixation Station for TB	\$44.00
SDL9750 .	IMMUNO-STAINING STATION SET	\$30.00

"5" on red bckgrnd=change
 PRICES EFFECTIVE 07/01/05

7/8/2005

2005 SDL Price List as of 7/1/05

SDL9752 .	IMMUNO STAINING RACK ONLY	\$22.50
SDL9755 .	REPLACEMENT CONTAINERS 10	\$35.00
SDL9780 .	PARAFFIN DISPENSER 1 Liter	\$250.00
SDL9785 .	PARAFFIN DISPENSER 2 Liter	\$300.00
SDL9801 .	BIOPSY SPONGE MATERIAL - 1000	\$37.00
SDL9803 .	MICROWAVE DRYING RACK-SMALL	\$32.50
SDL9804 .	Aqua Hold (each) replaces PAP PEN -SUPER	\$27.00
SDL9808 .	MICROWAVE DRYING RACK-LARGE	\$45.00
SDL9809 G	MICROWAVE JAR & RACK - GLASS	\$65.00
SDL9809 P	MICROWAVE JAR & RACK - PLASTIC	\$65.00
SDL9810 .	LABMAT MATERIAL/50 FT X 20"	\$40.00
SDL9812 .	BIO-HAZARD PLACEMAT (100)	\$110.00
SDL9815 .	BIOHAZARD WIPES (2000)	\$95.00
SDL9817 .	SUPERCHLOR (PACKAGE OF 100)	\$86.00
SDL9833	Soap for dispenser (12/case)	\$150.00
SDL9860 .	SLIDE RACK AND STAINING DISH - each	\$35.00
SDL9861 .	SLIDE RACK - each	\$10.00
SDL9862 .	STAINING DISH - each	\$30.00
SDL9863 .	SLIDE STORAGE BOX W/ 4 RACKS	\$40.00
SDL9865 .	LARGE SLIDE RACK - each	\$24.00
SDL9870 .	SPECIMEN CONTAINER W/CAP	\$41.81
SDL9871 .	STERILE CONTAINER W/CAP 1	\$14.69
SDL9930 .	SPLASH SHIELD - each	\$50.00
SDL9931 .	SPLASH GUARD each	\$30.00
SDL9932 .	SPLASH GUARD -each	\$30.00
SDL97100 .	ODOR CONTROLLER-TAN UNIT	\$42.50
SDL97101 .	ODOR CONTROLLER-SILVER DELUXE	\$45.00
SDL97106 .	ODOR CONTROLLER-SUPRA UNIT	\$56.00
SDL97110 .	ODOR CONTROLLER CARTRIDGE (5)	\$38.00
SDL97111 .	ODOR CONTROLLER CARTRIDGE (10)	\$73.50

Items highlighted in yellow and marked with * may be subject to minimum order.

All pricing is subject to change without notice. All prices are in US dollars and do not include freight.